

---

GSI 39

# GSI 39 Auto Tymp

---

## User Guide

January, 2009



## Safety summary

In this manual, two labels identify potentially dangerous or destructive conditions and procedures:

### **WARNING**

The **WARNING** label identifies conditions or practices that may present danger to the patient and/or user.

### **CAUTION**

The **CAUTION** label identifies conditions or practices that could result in damage to the equipment.

**NOTE:** Notes help you identify areas of possible confusion and avoid potential problems during system operation.

# **GSI 39**

---

## **Safety notes**

### **WARNING**

The GSI 39 is designed for compliance to IEC and UL 60601-1 when used in the patient vicinity. To achieve this compliance, use of hospital grade plug and receptacles is required. For patient and operator safety, the GSI 39 must be used with properly grounded plug and receptacles at all times. The GSI 39 is equipped with a specific power transformer (ref: 113-405800), which should not be interchanged with any other transformer or supply.


Any program aimed at obtaining reliable records of hearing thresholds should be staffed and supervised by appropriately-trained individuals.

Latex is not used anywhere in the manufacturing process.

The base material for the earphone cushions is made from natural and synthetic rubber.

The material used to manufacture Cardinal Health/GSI's eartips is Krayton Thermoplastic Rubber.

 **WARNING**

This symbol  indicates the location of a service adjustment part and is intended for service personnel only. The GSI 39 is a specifically calibrated audiometer and Tympanometer and, the periodic service and adjustments for the instrument that may be required should be done only by an authorized Cardinal Health/GSI service technician.

Please read the entire manual prior to using the GSI 39 to familiarize yourself with the test functions and proper accessory connections.

Accessory equipment connected to the analog and digital interfaces must be certified to the respective IEC standards (IEC950 for data processing or IEC 60601-1 for medical equipment). Furthermore, all configurations shall comply with the system standard IEC 60601-1-1. Everyone who connects additional equipment to the signal input or signal output port configures a medical system, and is therefore responsible that the system complies with the requirements of the system standard IEC60601-1-1. If in doubt, consult the technical service department or your local representative.

## Warranty

We, Cardinal Health/GSI warrant that this product is free from defects in material and workmanship and, when properly installed and used, will perform in accordance with applicable specifications. If within one year after original shipment it is found not to meet this standard, it will be repaired, or at our option, replaced at no charge except for transportation costs, when returned to an authorized Cardinal Health/GSI product service facility. If field service is requested, there will be no charge for labor or material; however, there will be a charge for travel expense at the service center's current rate.

**NOTE:** Changes in the product not approved in writing by Cardinal Health/GSI shall void this warranty. Cardinal Health/GSI shall not be liable for any indirect, special or consequential damages, even if notice has been given in advance of the possibility of such damages.

**THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OF FITNESS FOR A PARTICULAR PURPOSE.**

## Printing history

1st printing: November 2008

## **European authority representative**

Cardinal Health Germany 234 GmbH

Leibnizstrasse 7

97204 Hoechberg

Germany

Tel: (49) 931-4972 - 308

Fax: (49) 931-4972 - 62308

Email address: [support.nc.eu@cardinalhealth.com](mailto:support.nc.eu@cardinalhealth.com)



CE Mark per Medical Device Directive (93/42/EEC)

## Specifications

**Standards:** UL 60601-1 Medical Electrical Equipment Requirements for Safety  
IEC/EN 60601-1 General Requirements  
CSA C22.2 No.601-1-M90  
ANSI S3.39-1987 Aural Acoustic Impedance Admittance (Type 3)  
IEC 60645-5 Aural Acoustic Impedance/ Admittance (Type 3)  
ANSI S3.6-2004 Audiometers (Type 4)  
IEC 60645-1 Pure Tone Audiometers (Type 4)  
2004 Specifications for Audiometers (Type 4)  
PTB Certificate No. 15.11-94/53 Pure Tone Audiometers (Type 4)  
GL2005-00014 (ASHA 2005) Guidelines for Manual Pure-Tone  
Threshold Audiometry

**Protective Classification:** This system is intended for continuous operation and has a protective classification of Class 1.

**Class 1 Type B** 



---

## Tympanometry modes

<b>Probe Tone:</b>	<b>226 Hz, <math>\pm 2\%</math> 1000 Hz <math>\pm 2\%</math></b>
<b>Sound Pressure Level:</b>	<b>226 Hz:</b> 85.5 dB SPL, $\pm 2.0$ dB, measured in a 2.0 cm <sup>3</sup> coupler <b>1000 Hz:</b> 75 dB SPL, $\pm 2.0$ dB, measured in a 2.0 cm <sup>3</sup> coupler
<b>Harmonic Distortion:</b>	<3%
<b>Admittance (Compliance) Range:</b>	<b>226 Hz:</b> 0.0 to 1.5 cm <sup>3</sup> or 0.0 to 3.0 cm <sup>3</sup> <b>1000 Hz:</b> 0.0 to 5.0 mmho and 0.0 to 10 mmho

**NOTES:**

1. The range is automatically selected based upon the amplitude of the compensated tympanogram.
2. The maximum uncompensated (ECV + tympanogram peak) admittance (compliance) range is 0 to 5.0 cm<sup>3</sup>.
3. ECV/cavity limits for initiating pressurization is 0.2 to 5.0 cm<sup>3</sup>. Compliance Accuracy:  $\pm 0.1$  cm<sup>3</sup> or  $\pm 5\%$ , whichever is greater

## Pneumatic System

**Pressure Range:** +200 to -400 daPa

**NOTES:**

1. daPa = 1.02 mmH<sub>2</sub>O
2. For 226 Hz probe tone, pressure sweeps to at least -100 daPa. To save test time, pressure sweep stops once tympanogram returns to baseline after -100 daPa.
3. For 1000 Hz probe tone, pressure sweep does not stop until -400daPa.
4. Full pressure sweep for 5 cm<sup>3</sup> from sea level to 7000 ft. altitude with no leak.

**Pressure Accuracy:** ±10 daPa or ±15%, whichever is greater

**Rate of Sweep:** **226 Hz:** 600 daPa/sec except near tympanogram peak where sweep rate slows to 200 dapa/sec to provide better definition of peak compliance.

**1000 Hz:** 200 daPa/sec ±10 daPa/sec

**Direction of Sweep:** Positive to negative

**Tymp Test Time:** **226 Hz:** Approximately 1 second

**1000 Hz:** Approximately 3 seconds

**NOTE:** High compliance tympanograms will take somewhat longer.

**Gradient:** **226 Hz only:** Measurement of the tympanogram width taken at 50% of peak compliance.

## Acoustic Reflex Stimuli

<b>Frequencies:</b>	<b>226 Hz Probe Tone:</b> 500, 1000, 2000, and 4000 Hz for both ipsilateral and contralateral stimulation <b>1000 Hz Probe Tone:</b> 500, 2000 and 4000 Hz for both ipsilateral and contralateral stimulation
<b>Accuracy:</b>	$\pm 3\%$
<b>Total Harmonic Distortion:</b>	<5% for outputs less than 110 dBHL and <10% at 110 dBHL
<b>Rise/Fall Time:</b>	5 to 10 msec
<b>Transducers</b>	
<b>IPSI:</b>	Cardinal Health/GSI design
<b>CONTRA:</b>	Single Audiovox Model SM-N insert phone (Version 2 and 3 only)

## Output Levels:

### IPSI:

#### For 226 Hz Probe Tone:

**500 and 4000 Hz:** 80, 90, 100 db HL (Combo Probe design @ 4000 Hz has 80 and 90 db HL only)

**1000 and 2000 Hz:** 85, 95, 105 dB HL

#### For 1000 Hz Probe Tone:

**500 and 4000 Hz:** 80 and 90 dBHL

**2000 Hz:** 85 and 95 dBHL

### CONTRA:

#### For 226 Hz Probe Tone:

**500, 1000, 2000, 4000 Hz:** 90, 100, 110 db HL

#### For 1000 Hz Probe Tone:

**500, 2000 and 4000 Hz:** 90 and 100 dBHL

### NOTES:

1. **226 Hz Probe Tone:** Ipsi stimuli are time multiplexed with probe tone (93 ms ON, 66 ms OFF).  
**1000 Hz Probe Tone:** Ipsi stimuli are time multiplexed with probe tone (62 ms ON, 62 ms OFF).
2. Contra stimuli are steady tones.
3. Stimuli are presented at lowest level first. If there is no response, the intensity is increased by 10 dB until a response is detected or the maximum dB HL is reached.
4. Contra is available with Versions 2 and 3 only.

<b>Pressure:</b>	<b>226 Hz Probe Tone:</b> Reflex measures are set automatically to pressure at peak compliance with an offset of -20 daPa if peak pressure is negative and +20 daPa if peak pressure is positive. <b>1000 Hz Probe Tone:</b> Reflex measures are taken at 0 daPa regardless of peak pressure.
<b>Reflex Determination:</b>	<b>226 Hz:</b> Compliance change of 0.05 cm <sup>3</sup> or greater. <b>1000 Hz:</b> Compliance change of 0.1 mmho.
<b>Reflex Test Time:</b>	1 to 12 seconds depending upon the number of ipsi and/or contra test frequencies selected (4 maximum) and intensity required.

## Probe LED Indicators

### 226 Hz Probe:

<b>Steady yellow:</b>	Occlusion
<b>Blinking green:</b>	Ready to start
<b>Steady green:</b>	Test in progress
<b>Steady orange:</b>	Leak in pressure

### Combo Probe:

<b>Blinking green:</b>	Ready to start
<b>Steady green:</b>	Test in progress
<b>Steady orange:</b>	Occlusion
<b>Blinking orange:</b>	Leak in pressure
<b>No light:</b>	Test is finished

# GSI 39

---

## Audiometry mode (Versions 3 and 4 only)

**Frequencies:** 125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz

**Accuracy:**  $\pm 2\%$

**Total Harmonic Distortion:**  $< 2.5\%$  (125 to 3000 Hz measured acoustically at maximum dB HL; 4000 and 6000 Hz measured electrically)

## Transducers

**Audiometric Headset:** Pair TDH-39 earphones with MX41AR cushions (60 ohms impedance) - Versions 3 & 4 only, Headband force per ANSI S3.6 and IEC 645 ( $4.5 \pm 0.5$ )

**Insert Earphones:** ER-3A or ER-5A Insert earphones (50 ohm impedance)

---

## Intensity Levels

**TDH-39 Headphones**

<b>125 Hz</b>	-10 to 50 dB HL
<b>500 to 6000 Hz</b>	-10 to 90 dB HL
<b>250 and 8000 Hz</b>	-10 to 70 dB HL

**Insert Phones**

<b>125 Hz</b>	-10 to 40 dBHL
<b>500 to 4000 Hz</b>	-10 to 80 dBHL
<b>6000 Hz</b>	-10 to 70 dBHL
<b>250 and 8000 Hz</b>	-10 to 60 dBHL

**NOTE:** An additional +10 dB is available per frequency via the +10 dB button.

<b>Accuracy:</b>	125 to 4000 Hz $\pm 3$ dB 6000 and 8000 Hz $\pm 5$ dB
<b>Step Size:</b>	5 dB
<b>Signal-to-Noise Ratio:</b>	> 70 dB in 1/3 octave; less than -10 dB HL for levels less than 60 dB HL
<b>Rise/Fall Time:</b>	20 to 50 msec

## Tone Format

Tone is normally off until the Present bar is depressed.

<b>Continuous</b>	Tone is steady when present bar is depressed
<b>Pulsed</b>	Tone is pulsed at 2.5/sec (i.e., 200 msec ON, 200 msec OFF)
<b>FM (frequency modulated)</b>	Tone is frequency modulated at a rate of 5 Hz, $\pm 5\%$

# GSI 39

---

## Printer

<b>Paper Roll Length:</b>	Approximately 80 feet (960")
<b>Tests/Roll:</b>	
<b>Versions 1 and 2:</b>	Approximately 420 Tymps/Reflex or 210 people
<b>Versions 3 and 4:</b>	Approximately 230 tests or 115 people
<b>Assumption:</b>	2 Tymps/Reflex + 1 Audiogram per person
<b>Speed:</b>	Approximately 1 minute to print three screens: Tympanogram Tympanogram + reflex (4) Audiogram
<b>External Printer:</b>	Optional Deskjet color printer recognizing PCL3 or PCL3 GUI; 8-1/2" x 11" or A4 format

## Power

<b>Line Voltage:</b>	100 - 240 VAC ( $\pm 10\%$ ) <b>NOTE:</b> Desktop power supply.
<b>Frequency Range:</b>	47 - 63 Hz ( $\pm 5\%$ )
<b>Power Consumption:</b>	16 watts maximum while printing. Low voltage input for desktop power supplies 7 VDC, 5.0 A.
<b>Display:</b>	240 x 64 graphical, monochrome LCD



## Environmental

<b>Temperature:</b>	
<b>Operating:</b>	59° to 104° F (15° to 40° C)
<b>Warm-up time:</b>	10 minutes for instruments stored at room temperature.
<b>Storage/Shipping:</b>	-93° to 149° F (-69° to 65° C)
<b>Ambient Pressure:</b>	98 kPa to 104 kPa
<b>Humidity:</b>	15% to 95%

## Mechanical - Instrument

### Instrument

<b>Dimensions:</b>	12.5" W x 14.5" D x 4.7" H 31.8 cm W x 36.8 cm D x 11.9 cm H
<b>Weight:</b>	5 lbs (2.3 kg) - unit and probe

### Shipping Carton

<b>Dimensions:</b>	19.5" W x 22.5" D x 8.25" H 49.5 cm W x 57.2 cm D x 20.9 cm H
<b>Weight:</b>	13.1 lbs (6 kg)

# **GSI 39**

---

## **Supplied Accessories**

### **1739-9710 GSI 39**

Probe (226 Hz version)	1739-3200
Contra Insert Phone (Versions 2 and 3 only)	8000-0079
TDH-39 Headset (Versions 3 and 4 only)	8000-0178
Test Cavity (226 Hz)	1700-1030
Power Supply	113-405800
Eartips, (Probe) 6 sizes, 2 each (all versions)	1700-9622
Eartips, (Contra Insert Phone) color coded, 8 sizes, 4 each (Versions 2 & 3 only)	1700-9660
Paper 3 rolls, thermal 4" wide (all versions). Shipped with system.	
Note: To reorder a single roll of paper, use the following part number:	1738-9601
Instruction Manual (all versions)	1739-0100
GSI 39 Quick Reference Guide	1739-0160

### **1739-9711 GSI 39**

Combo Probe (226 Hz / 1 KHz version)	1739-3250
Contra Insert Phone (Versions 2 and 3 only)	8000-0079
TDH-39 Headset (Versions 3 and 4 only)	8000-0178
Test Cavity (226 Hz)	1700-1030
Power Supply	113-405800
Eartips, (Probe) 6 sizes, 2 each (all versions)	1700-9622
Eartips, (Contra Insert Phone) color coded, 8 sizes, 4 each (Versions 2 & 3 only)	1700-9660
Paper 3 rolls, thermal 4" (10.16 cm) wide (all versions). Shipped with system.	
Note: To reorder a single roll of paper, use the following part number:	1738-9601
Instruction Manual (all versions)	1739-0100
GSI 39 Quick Reference Guide	1739-0160
Probe Cleaning Kit	2000-9610
Probe Mount - Shoulder	1700-9646
Probe Mount - Wrist	1700-9642
Probe Mount - Clothes	1700-9608

## Optional Accessories

Dust Cover	1738-9620
Carrying Case	1738-9680
Patch Cord (1)	4204-0505
Subject Response Handswitch	7874-0156
Audiocups Earphone Sound Enclosures	8000-0155
Service Manual	1739-0110
Insert Phone Assembly 3A (50 ohm impedance)	1700-0441
Insert Phone Assembly 5A (50 ohm impedance)	1700-0882



### **WARNING**

Use only Cardinal Health/GSI supplied components and accessories.

## Catalog Listings: GSI 39 Auto Tympanometer

		Tymp + Ipsi V.1	Tymp + Ipsi + Contra V.2	Tymp + Ipsi + Contra + Audio V.3	Tymp + Ipsi + Audio V.4	Tymp only V.5
226 Hz only	USA	1739-9710	1739-9720	1739-9730	1739-9740	1739-9750
	International	1739-9710INT	1739-9720INT	1739-9730INT	1739-9740INT	1739-9750INT
Combo 226 + 1k Hz	USA	1739-9711	1739-9721	1739-9731	1739-9741	1739-9751
	International	1739-9711INT	1739-9721INT	1739-9731INT	1739-9741INT	1739-9751INT

## **Country Kits**

(one must accompany each system)

**Each Country Kit contains a Power Cord, User Manual, Wall Chart and a Quick Reference Guide**

1739-9400	North America cord, English
1739-9401	CtryKit, North America Cord., Spanish
1739-9402	CtryKit, North America Cord, Portuguese
1739-9403	CtryKit, North America Cord., French
1739-9404	CtryKit, Europe Cord, French
1739-9405	CtryKit, Europe Cord, German
1739-9406	CtryKit, Europe Cord, Spanish
1739-9407	CtryKit, Europe Cord, Portuguese
1739-9408	CtryKit, Europe Cord, Russian
1739-9409	CtryKit, Europe Cord, English
1739-9410	CtryKit, UK-Ireland Cord, English
1739-9411	CtryKit, Italy Cord, Italian
1739-9412	CtryKit, Italy Cord, Spanish
1739-9413	CtryKit, Swiss Cord, German
1739-9414	CtryKit, Swiss Cord, French
1739-9415	CtryKit, Swiss Cord, English
1739-9416	CtryKit, Denmark Cord, English
1739-9417	CtryKit, Israel Cord, English
1739-9418	CtryKit, South Africa-India Cord, English
1739-9419	CtryKit, Australia Cord, English
1739-9420	CtryKit, Australia Cord, Chinese
1739-9421	CtryKit, North America Europe Cord, Korean
1739-9422	CtryKit, North America Cord, Japanese

## Upgrade Kits

1739-9600 \*Upgrade from 226 Hz to 226 Hz + 1000 Hz probe (use same number for all versions)

1739-9612 \*Upgrade GSI 39: 226 Hz Version 1 to 226 Hz Version 2

1739-9613 \*Upgrade GSI 39: 226 Hz Version 1 to 226 Hz Version 3

1739-9614 \*Upgrade GSI 39: 226 Hz Version 1 to 226 Hz Version 4

1739-9623 \*Upgrade GSI 39: 226 Hz Version 2 to 226 Hz Version 3

1739-9643 \*Upgrade GSI 39: 226 Hz Version 4 to 226 Hz Version 3

1739-9651 \*Upgrade GSI 39: 226 Hz Version 5 to 226 Hz Version 1

1739-9652 \*Upgrade GSI 39: 226 Hz Version 5 to 226 Hz Version 2

1739-9653 \*Upgrade GSI 39: 226 Hz Version 5 to 226 Hz Version 3

1739-9654 \*Upgrade GSI 39: 226 Hz Version 5 to 226 Hz Version 4

(\* All upgrade kits require calibration. Consult your local Cardinal Health/GSI Distributor for pricing.)

Blank page.

<b>Safety summary</b> .....	<b>a</b>
Safety notes.....	b
<b>Warranty</b> .....	<b>d</b>
<b>Printing history</b> .....	<b>d</b>
<b>European authority representative</b> .....	<b>e</b>
<b>Specifications</b> .....	<b>f</b>
Tympanometry modes.....	g
Pneumatic System.....	h
Acoustic Reflex Stimuli.....	i
Probe LED Indicators .....	k
Audiometry mode (Versions 3 and 4 only).....	l
Transducers .....	l
Intensity Levels.....	m
Tone Format .....	m
Printer.....	n
Power .....	n
Environmental.....	o
Mechanical - Instrument .....	o
Supplied Accessories .....	p
Optional Accessories .....	q
Catalog Listings: GSI 39 Auto Tymp .....	q
Country Kits.....	r
Upgrade Kits .....	s

## **Introduction 1**

<b>Introduction</b> .....	<b>1-3</b>
<b>Customer responsibility</b> .....	<b>1-4</b>
Tympanometry and Gradient .....	1-6
Gradient .....	1-8
Screening acoustic reflex .....	1-9
Screening audiometry .....	1-11
<b>Recycling / disposal</b> .....	<b>1-12</b>
<b>Glossary of terms</b> .....	<b>1-13</b>

## **Installation 2**

<b>Unpacking and Inspection</b> .....	2-3
<b>Printer and Display</b> .....	2-4
<b>Rear Panel Labels and Connectors</b> .....	2-5
<b>Bottom Panel</b> .....	2-6
Table of symbols on the GSI 39.....	2-7
<b>Initial set-up</b> .....	2-8
Loading the paper .....	2-9
Paper storage.....	2-10

## **Operation 3**

<b>226 Hz Probe Indicators</b> .....	3-3
<b>Combo Probe Indicators (226 Hz and 1000 Hz probe tone)</b> .....	3-4
Preparing the probe assembly .....	3-4
<b>Front Panel Controls and Indicators</b> .....	3-6
<b>Tympanometry testing information</b> .....	3-13
Helpful hints .....	3-13
Obtaining a seal .....	3-14
<b>Combo Probe Insertion</b> .....	3-17
<b>Audiometry testing information (Versions 3 and 4)</b> .....	3-18
Instructing the patient/subject.....	3-18
Placement of earphones .....	3-19
Placement of Insert earphones .....	3-19
Response handswitch (optional accessory).....	3-19
<b>Tympanometry/Reflex Test Sequence</b> .....	3-20
a. Tympanometry only mode .....	3-20
b. Tympanometry and Ipsilateral Reflex.....	3-23
c. Temporary programming of ipsilateral acoustic reflex test frequencies .....	3-26
d. Tympanometry and Contralateral Reflex (Version 2 and 3).....	3-27
e. Tympanometry and ipsilateral/contralateral reflexes (versions 2 and 3 only) .....	3-29
Exit tympanometry/reflex .....	3-30
<b>Audiometry Test Sequence (Versions 3 and 4 only)</b> .....	3-31
To enter the Audiometry mode .....	3-31
Transducer Selection.....	3-31
To change the frequency .....	3-33
To change the intensity level of the test tone.....	3-34
Screening audiometry .....	3-36
<b>Audiometric Threshold</b> .....	3-37
Manual Threshold Audiometry.....	3-37
Automatic Hearing Level.....	3-39
Performing the Auto HL Procedure .....	3-40



---

Exit audiometry.....	3-40
<b>Tests in memory .....</b>	<b>3-41</b>
Memory erase .....	3-41
Printing test results.....	3-42

## **Program Mode 4**

<b>Program Mode .....</b>	<b>4-3</b>
Program Mode Menu Items .....	4-4
<b>Program Menu Page 1 Option Descriptions.....</b>	<b>4-5</b>
PROBE HZ . . . . .	4-5
TYMP OPTIONS . . . . .	4-5
REFLEX DISPLAY .....	4-8
• Reflex dB HL plus curve.....	4-8
• Reflex dB HL only .....	4-8
• Reflex yes/no.....	4-9
226 Hz REFLEX.....	4-10
1000 Hz REFLEX.....	4-10
AUTO HL SETUP.....	4-11
LANGUAGE .....	4-14
AUD RANGE NORMAL/AUD RANGE NARROW .....	4-14
PRINT – AUDIOGRAM / PRINT – AUD TABLE .....	4-15
DEF XDUCER TDH39 / DEF XDUCER INSERT .....	4-15
<b>Program Menu Page 2 option descriptions .....</b>	<b>4-16</b>
DATA XFER CONFIG .....	4-16
POWER UP SETTINGS .....	4-16
PRN HEADER GSI/PRN HEADER OFF/PRN HEADER CUSTOM.....	4-17
INTERNAL PRINTER/EXTERNAL PRINTER.....	4-18
RESET TO DEFAULTS .....	4-19
Exiting the program mode .....	4-20

**Routine Maintenance 5**

**PreTest Tymp checks ..... 5-3**  
Calibration Quick Check for 226 Hz Probe ..... 5-4  
Calibration Quick Check for Combo Probe ..... 5-5  
Altitude adjustment ..... 5-7  
**Pre-Test Audiometric Checks (Version 3 and 4 only) ..... 5-9**  
Noise recovery period ..... 5-9  
Elimination of ambient noise ..... 5-10  
**Biological Check ..... 5-11**  
**Preventive Maintenance ..... 5-12**  
Cleaning the system ..... 5-12  
Recommended cleaning solutions ..... 5-12  
Cleaning patient contact reusable devices ..... 5-13  
**Probe care - 226 Hz Probe ..... 5-14**  
Probe nose cone cleaning ..... 5-14  
**Probe Care - Combo Probe Tip ..... 5-17**  
**Earphone Care (Versions 3 and 4 only) ..... 5-19**  
**Paper supply ..... 5-20**

**Test Results 6**

**Ear Canal Volume - 226Hz Probe Tone ..... 6-3**  
Normal ..... 6-3  
Abnormal ..... 6-3  
**Compliance Peak ..... 6-4**  
Normal ..... 6-4  
Abnormal ..... 6-4  
**Pressure Peak ..... 6-5**  
Normal ..... 6-5  
Abnormal ..... 6-5  
**Gradient ..... 6-5**  
**Acoustic reflex ..... 6-6**  
Normal ..... 6-6  
Abnormal ..... 6-6  
**Audiometry ..... 6-6**  
**Special Messages and Error Codes ..... 6-7**  
**Sample Test Results ..... 6-8**

**Computer Interface 7**

**Introduction**..... 7-3  
**Operation**..... 7-3  
Other LCD screen messages ..... 7-4  
**Data Transfer Program Mode** ..... 7-5  
**Computer Interface** ..... 7-6  
Interface configuration..... 7-6  
Cable connections ..... 7-6

**Electromagnetic Compatibility (EMC) 8**

GSI 39 - Electromagnetic Compatibility (EMC) Information..... 8-3

**Bibliography 9**

Blank page.

# Chapter 1

# Introduction

---

Blank page.

## Introduction

The **GSI 39 Auto Tymp** (hereafter referred to as '**instrument**' in this guide unless otherwise noted for clarity) is a versatile combination instrument that provides testing capability for tympanometry alone, tympanometry combined with screening acoustic reflex measurements, and screening audiometry.

Five different versions are available to meet your individual testing needs.

- Version 1 provides two modes of operation, tympanometry alone and tympanometry plus screening ipsilateral acoustic reflex testing.
- Version 2 permits tympanometry alone and tympanometry combined with ipsilateral and contralateral screening acoustic reflex measurements.
- Version 3 provides testing capability for all three test modes i.e., tympanometry alone, tympanometry combined with ipsilateral and contralateral screening acoustic reflex measurements, and screening audiometry, both manual and automated.
- Version 4 allows tympanometry alone, tympanometry combined with ipsilateral acoustic reflex screening testing, and manual and automated screening audiometry.
- Version 5 is for tympanometry only.

It is possible to field retrofit versions 1, 2, 4 and 5 with the full functionality provided with version 3 after the time of original purchase.

Each version can be fitted with the Combo Probe that allows for both 226 Hz and 1000 Hz probe tone. 1000 Hz probe tone is recommended for testing on infants 0 - 6 months of age.

An optional soft-sided carrying case is available for portability. Also, a patient handswitch, patch cords, and earphone sound enclosures may be purchased as optional accessories.

## Customer responsibility

### **WARNING**

This product and its components will perform reliably only when operated and maintained in accordance with the instructions contained in this manual, accompanying labels, and/or inserts. A defective product should not be used. Make sure all connections to external accessories are snug and secured properly. Parts which may be broken or missing or are plainly worn, distorted or contaminated should be replaced immediately with clean, genuine replacement parts manufactured by or available from Cardinal Health/GSI.

This product should not be used in the presence of fluid that can come into contact with any of the electronic components or wiring. Should the user suspect fluids have contacted the system components or accessories, the unit should not be used until deemed safe by a Cardinal Health/GSI certified service technician.

Do NOT use in the presence of flammable gaseous mixtures. Users should consider the possibility of explosions or fire when using this device in close proximity to flammable anesthetic gases.

Periodically have a service technician perform electrical safety checks on the unit in order to show continued compliance to IEC and UL 60601-1.





**WARNING** The GSI 39 is designed to be used with a hospital grade outlet. Injury to personnel or damage to equipment can result when a three-prong to two-prong adapter is connected between the GSI 39 power plug and an AC outlet or extension cord. Additionally, the GSI 39 is equipped with a specific power transformer (113-405800), which should not be interchanged with any other transformer or supply.

The GSI 39 is a specifically calibrated device and the periodic service and adjustments, which the instrument may require, should be done only by an authorized Cardinal Health/GSI service technician.



**CAUTION**

The GSI 39 is designed to comply with the EMC requirements according to IEC 60601-1-2. Radio transmitting equipment, cellular phones, etc. shall not be used in the close proximity of the device since this could influence the performance of the device. Particular precaution must be considered during use of strong emission sources such as High Frequency surgical equipment and similar device so that, for example, the HF-cables are not routed on or near the device. If in doubt, contact a qualified technician or your local representative.

## Tympanometry and Gradient

Tympanometry is an objective technique used since the late 1960's to measure the mobility (compliance) and the pressure within the middle-ear system. During tympanometry, a tone (e.g., 226 Hz probe tone) is presented to the ear canal via the light-weight probe. The probe tone is used to measure the compliance changes within the middle-ear system while air pressure within the hermetically sealed ear canal is varied from a positive to a negative value. A positive pressure within the ear canal space, in the absence of middle-ear pathology, causes the middle-ear system to stiffen up or become less mobile. This is caused by the pressure difference between the sealed ear canal space and the middle-ear space, which forces the tympanic membrane to stretch inward. A stiffened middle-ear system displays little or no compliance. As the pressure within the ear canal is brought back toward atmospheric (ambient or 0 daPa) pressure, the pressure difference between the ear canal space and the middle-ear space is reduced in normal ears.

At or near atmospheric pressure (0 daPa), the greatest amount of sound (probe tone) enters the middle-ear system. This is the air pressure value where the middle-ear system displays the maximum amount of compliance.

When the air pressure within the ear canal is then reduced to a negative value with respect to the middle-ear space, a pressure difference is once again established and the middle-ear system becomes less compliant. Therefore, by varying the pressure within the ear canal, it is possible to make a series of compliance measurements by means of the probe tone. The tracing which depicts these compliance changes is referred to as a tympanogram. The point of the tympanogram which represents the point of maximum compliance is the compliance peak of the tympanogram. The air pressure (pressure at the peak) where this compliance peak occurs approximates the pressure within the middle-ear system, since maximum mobility is only possible when there is little or no pressure difference between the ear canal and the middle-ear space. Compliance using a 226 Hz probe tone is measured with respect to the ability of an equivalent volume of air to conduct sound and the scientific quantity used is  $\text{cm}^3$ . Compliance using the 1000 Hz probe tone is measured in mmho's as it is not an ear canal compensated measurement. Air pressure is measured in decaPascals (daPa).

**NOTE:**  $1.02 \text{ mm H}_2\text{O} = 1.0 \text{ daPa}$

The presence of a pathological condition which interferes with the mobility of the tympanic membrane, the ossicular chain, or the air pressure within the middle-ear space can be detected during tympanometry. For example:

- If the air pressure within the middle-ear space becomes negative due to a blocked eustachian tube, tympanometry permits the measurement of this negative pressure and its effect on middle-ear compliance.
- If fluid builds up within the middle-ear space, this fluid will restrict the ability of the ossicular chain to conduct sound to the cochlea. If small air pockets exist within the fluid, the tympanogram will indicate the negative pressure where the restricted mobility occurs. With a totally fluid-filled middle-ear space, no mobility will be measured during tympanometry at any pressure value.
- In the case of a “glue-ear”, the ossicular chain is restricted in mobility. This tympanogram would depict a flat line with no identifiable pressure peak.

## Gradient

Gradient (width) measurements are used to describe the shape of a tympanogram in the vicinity of the peak. Often, the presence or absence of fluid in the middle ear is not clearly indicated by otoscopy and tympanometry alone. This evaluation is especially difficult when the peak pressure is in the normal range.

The presence of fluid within the middle-ear space alters the shape of a tympanogram (i.e., makes the tympanogram wider near its peak). A larger-than-normal gradient can indicate the presence of fluid in the middle ear when other parameters are within normal limits. In this way, the gradient acts as an adjunct to the tymp and ear canal volume measurements by helping to differentiate between tymps with similar peak values.

The instrument uses tympanometric width to determine the gradient by measuring the pressure interval at one-half of the tymp peak height. Differing tymp peak widths can point to different middle-ear conditions, even when peak height and pressure are within normal range. For example, middle-ear effusion brought on by secretory otitis media might result in an increased tympanogram width and, therefore, an increased gradient value. This would occur because the ossicular chain cannot react to the change in pressure introduced during the tympanogram in the same way that it would if the middle ear were properly aerated. The continued presence of effusion, leading eventually to a completely fluid filled middle-ear cavity, will reduce the magnitude of the tympanogram to the point where no change in compliance is detectable across the pressure range. Under this condition, no gradient measurement is possible.

On the GSI 39, gradient measures are only calculated for the 226 Hz probe tone conditions.

## Screening acoustic reflex

An acoustic reflex occurs when a very loud sound (stimulus) is presented to the auditory pathway. During acoustic reflex testing, the stimulus is presented to the ear canal through a probe (ipsilateral) or through an insert phone (contralateral). This stimulus then travels through the middle ear to the cochlea. From the cochlea, frequency and intensity information is transmitted via the 8th nerve to the brain stem where a determination is made as to whether or not the intensity of the stimulus is high enough to elicit the reflex response. If it is, a bilateral response occurs (i.e., the right and left 7th nerves innervate their respective middle-ear muscles (stapedial muscles) causing them to contract). As these muscles contract, they stiffen their respective ossicular chains. This stiffening of the ossicular chain reduces the compliance of each middle-ear system. As in tympanometry, a probe tone is used to measure this decrease in compliance.

When the stimulus is presented to the same ear where the measurement takes place, the test is referred to as an ipsilateral (same side) acoustic reflex test. When the stimulus is presented to the opposite ear from where the measurement takes place, the test is referred to as a contralateral (opposite) acoustic reflex test.

During ipsilateral acoustic reflex testing, both the stimulus and the probe tone are presented via the hand-held probe. With contralateral testing, the stimulus is presented via an insert phone or earphone and the probe tone is presented via the hand-held probe. In both cases, the measurement is made from the ear where the probe is positioned. For 226 Hz probe tone reflex measurements, the air pressure within the ear canal where the probe is positioned is set to the pressure value measured at the point of maximum compliance for that ear during tympanometry with an offset of -20 daPa (or +20 daPa for a positive pressure peak).

For 1000 Hz probe tone reflex measurements, the system will measure the change in compliance at 0 daPa, regardless of peak pressure. When evaluating the absence of reflexes, you should note the peak pressure of the tympanogram. Decreased mobility at 0 daPa may explain the lack of reflexes. Reflex testing should be repeated when the middle ear pressure has returned to 0 daPa.

Acoustic reflex measurements are useful to determine the integrity of the neuronal pathway involving the 8th nerve, brainstem, and the 7th nerve. Since the acoustic reflex test (ipsilateral or contralateral) is performed at high intensity levels and since it involves a measurement of middle-ear mobility, acoustic reflex testing is not a test of hearing.

The acoustic reflex also serves as a good validation of tympanometric results since an acoustic reflex cannot be measured in the absence of a compliance peak. In other words, if the tympanometric results indicate no mobility over the pressure range available with your instrument, no reflex can be measured. If the test results indicate a reflex response in the absence of a compliance peak, one has cause to question the validity of the tympanometric test results. This indicates that the tympanogram should be repeated.

Clinical middle-ear instruments allow the measurement of the acoustic reflex threshold since they provide the ability to manually change the intensity of the stimulus to a level where a reflex response is just barely detectable for each patient tested. However, the instrument automatically presents the stimulus in a very definite stimulus intensity sequence. This preset intensity sequence may start at a level above an individual's acoustic reflex threshold level. Also, since the instrument uses a hand-held probe and noise from hand motion can be detected by the instruments circuitry, the magnitude of a detectable response must be somewhat higher than the criterion generally used during clinical acoustic reflex threshold testing to avoid artifact caused by hand motion. Thus, the acoustic reflex measurements made with the instrument are referred to as screening acoustic reflex testing. The purpose of these screening reflex tests is to determine whether a reflex is detectable rather than to determine the lowest intensity at which the reflex occurs (i.e., threshold testing).

## Screening audiometry

While tympanometry and acoustic reflex measurements check the integrity of the middle-ear system, audiometry provides a means for checking the integrity of the entire auditory pathway. Screening audiometry provides a method to check an individual's ability to hear a test signal at a particular intensity level or at the lowest possible intensity level without the use of masking.

During screening audiometry, the test signal is generally presented through an earphone to the ear under test. Different screening test protocols define the frequencies and intensity sequence to be used to obtain a response. Audiometric testing requires a behavioral response for the individual being tested. This consists of having the individual raise a finger/hand or press a handswitch (optional) whenever the test signal is heard. The finger/hand is lowered or the handswitch is released when the test signal is no longer audible. Thus, the individual being tested must be able to understand a set of simple instructions and have the ability to provide some physical sign when the test signal is heard.

The GSI 39 allows for both manual and automated audiometry. For further details on automated audiometry, see *Automatic Hearing Level* in Chapter 3 of this guide.

## Recycling / disposal



**Many local laws and regulations consider electric equipment-related waste as hazardous or requiring special procedures for recycling or disposal. This includes batteries, printed circuit boards, electronic components, wiring and other elements of electronic devices. Follow all of your respective local laws and regulations for the proper disposal of batteries and any other parts of your system.**

Check the Cardinal Health website for recommended instructions and addresses for proper return or disposal of electronic wastes relating to Cardinal Health/GSI products in Europe and other localities.

The contact information for the WEEE - In Europe

Cardinal Health GmbH  
NeuroCare Group  
D-97204 Hoechberg  
Germany

Tel: +49 (0) 931-4972-309

Fax: +49 (0) 931-4972-318



## Glossary of terms

**Acoustic Reflex** - Reflex arc elicited in the presence of very loud sounds, which cause a decrease in middle-ear compliance as a protective mechanism for the cochlea.

**Automated Audiometry (Auto HL)** - The Audiometer uses an adaptive computerized procedure that allows the listener to control the intensity with a hand switch.

**Compliance Peak** - The point of maximum mobility in a tympanogram, which indicates the degree of mobility within the middle-ear system.

**Contralateral Acoustic Reflex** - The acoustic reflex elicited when the stimulus is presented to the opposite ear from where the response is measured.

**Ear Canal Volume** - Volume measured between the tip of the probe and the tympanic membrane at the starting pressure for a tympanogram using a 226 Hz probe tone.

**Ipsilateral Acoustic Reflex** - The acoustic reflex elicited when the stimulus is presented to the same ear where the response is measured.

**Normal Box** - Range of pressure peak and compliance peak values associated with normal middle-ear function. (-150 daPa to +100 daPa, 0.2 cm<sup>3</sup> to 1.4 cm<sup>3</sup> per ASHA, 32, Supl. 2, 1990, 17-24) - only available on 226 Hz probe tone testing.

**Pressure Peak** - Pressure value where maximum mobility occurs in a tympanogram. This pressure value approximates the pressure within the middle-ear space.

**Probe Tone** - Tone used to measure middle-ear mobility (226 Hz or 1 kHz).

**Screening Audiometry** - A hearing test performed with a variety of frequencies and intensities without the use of masking to determine if an individual can hear.

**Tympanogram** - The tracing which depicts the results of tympanometry.

**Tympanometry** - An objective measurement of middle-ear mobility and middle-ear pressure through the use of a low frequency sound (probe tone) and air pressure changes.

Blank page.

# Chapter 2

# Installation

---

Blank page.

## Unpacking and Inspection

Examine the outside of the shipping container for any signs of damage. Notify your carrier immediately if any damage is noted.

Carefully remove your instrument from its shipping container. Remove the plastic bag protecting the instrument. If the instrument appears to have suffered mechanical damage, notify the carrier immediately so that a proper claim can be made. Be certain to save all packing material so that the claim adjuster can inspect it as well. As soon as the carrier has completed the inspection, notify your Cardinal Health/GSI Distributor.

Check that all accessories listed in Table 1 (per version ordered) are received in good condition. If any accessories are missing or damaged, notify your Cardinal Health/GSI Distributor or the factory immediately. See the *Specifications* listed in the first section of this manual for the catalog numbers of accessories and also for a listing of optional accessories.

<b>Probe Assembly</b>	<b>Instruction Manual</b>
<b>Power Module</b>	<b>Contralateral Insert Phone (Versions 2 and 3)</b>
<b>Test Headset (Versions 3 and 4 only)</b>	<b>Probe Eartips (6 sizes, 2 each)</b>
<b>Eartips for Insert Phone (8 sizes, 4 each) (Versions 2 and 3 only)</b>	<b>Paper (3 rolls)</b>
<b>Test Cavity</b>	

**Table 1: Accessories supplied**

**NOTE:**Keep the original packing material and shipping container so the instrument can be well packaged if it needs to be returned to the local service center for repair or calibration.




### **WARNING**


**Only use Cardinal Health/GSI approved parts and accessories.**

**NOTE:**The use of parts or materials that are not recognized to be used with the device can degrade minimum safety.


**Printer and Display**



**Figure 1: GSI 39 major component identifications.**



**Figure 2: Probe (226 Hz only).**



**Figure 3: Combo Probe (226 Hz / 1 KHz)**

## Rear Panel Labels and Connectors





Figure 4: Rear panel labels and connectors.

- R1:** Symbol denotes a Type B.
- R2:** Symbol denotes Attention, consult accompanying documents.
- R3:** Connector for handswitch.
- R4:** Connector for contralateral insert phone.
- R5:** Connectors for right and left earphone.
- R6:** Power Input Jack for external power supply.
- R7:** Power Switch with ON/OFF indicators.
- R8:** USB port for connecting to an external printer.
- R9:** USB port for connecting to a computer.

**! WARNING** Accessory equipment connected to the analog and digital interfaces must be certified to the respective IEC standards (IEC 950 for data processing or IEC 60601-1 for medical equipment). Furthermore, all configurations shall comply with the system standard IEC 60601-1-1. Everyone who connects additional equipment to the signal input or signal output part configures a medical system, and is, therefore, responsible that the system complies with the requirements of the system standard IEC 60601-1-1. If in doubt, consult the technical service department or your local representative.















**Bottom Panel**



**Figure 5: Bottom panel.**



**Table of symbols on the GSI 39**

No.	Symbol	Description
1		Attention, consult accompanying documents
2		Date of manufacture
3		CE Marked in accordance with the European Council Directive 93/42/EEC concerning medical devices
4		Medical device listing mark for U.S. and Canada by Intertek Testing Service
5		Special Recycling Required. Do not dispose in landfill
6		Type B equipment
7	REF	Symbol for “CATALOG NUMBER”
8		Stand-by
9		Right Ear
10		Left Ear
11		Patient Response Button
12		AC Power
13		Printer Connector
14		USB Type Connectors
15		Computer Connector

## Initial set-up

Place the instrument on a stable counter or table where it will subsequently be used. The location should be near a properly grounded wall outlet. Carefully attach purchased accessories to their appropriately labelled connector on the rear panel of the instrument (see Figure 4 on page 2-5).

Locate the **power** switch on the rear panel of the instrument and move the switch to the **On** position. Once power is turned on, the light on the LCD will be illuminated and the orange light on the probe will be lit. The display on the LCD will display a scroll bar across the top to indicate the system is initializing.

The system will power up to the factory default test mode (to set user-defined power up setting, see *Program Mode* in Chapter 4) and the probe green lamp will begin to blink indicating that the instrument is ready to begin the tympanometry. If both the green and yellow lamps are illuminated at the same time following power on, the probe is occluded or the tympanometry software did not initialize properly. Simply move the power switch to the off position, inspect the probe tip for any signs of an occlusion, and reposition the power switch to **On**. If both green and yellow lamps are still illuminated and you are certain that the probe is not occluded, contact your local service representative or the Cardinal Health/GSI service department for repair. In the mean time, it is still possible to use the Audiometry mode (if purchased).

Allow the instrument to warm-up for about 10 minutes before conducting a test. This allows the electronic circuits to stabilize prior to use. If the storage temperature is lower than the room temperature, allow some additional time for the instrument to reach room temperature.




### **CAUTION**


**Use only the Cardinal Health/GSI provided power supply. The GSI 39 provided power supply should only be connected to a power source meeting the following range 90-246VAC, 47-63Hz. In North America the power source should be a maximum of 120VAC.**

## Loading the paper

Remove the printer cover by placing your fingers along the back edge of the printer and pulling upward on the cover. Cut the printer paper so that the leading edge of paper is straight across. Place the roll of paper inside the paper well so that the paper will unroll from the lower surface. See the paper loading label located on the side of the paper well.



**Figure 6: Paper loading.**

Position the leading edge of the paper roll into the paper slot. Press the paper advance  button until a section of paper is long enough to pass through the printer cover.

## **Paper storage**

The instrument is supplied with a thermal printer. This type of printer requires a heat-sensitive paper to create an image. For maximum paper life, any spare rolls of paper should be stored as follows:

- a. Store in the dark (i.e., in a drawer or cabinet)
- b. Do not store above 77° F (25° C)
- c. Store at less than 65% relative humidity

The above recommendations are for the maximum paper life (greater than five years). Storing your thermal paper at high temperatures or high humidity levels will shorten the total paper life, depending on the actual temperature and humidity to which the paper is subjected. The paper will show some darkening if stored for 24 hours at 113° F (45° degrees C) and a relative humidity of greater than 90%, so avoid leaving your paper in a hot car or other hot area overnight. Always avoid storing unused paper or printed tests in a lighted area.


# Chapter 3 Operation

---

Blank page.

---

## 226 Hz Probe Indicators




**Figure 1: Probe indicators**

- P1 - Yellow:** The probe is occluded. Remove the probe and inspect for cause of occlusion.
- P2 - Green lamp:** *Blinking* - The instrument is ready to begin a Tymp.  
*Steady green* - Test successfully started and in progress.
- P3 - Orange:** A pressure leak has been detected.

## Combo Probe Indicators (226 Hz and 1000 Hz probe tone)

### Preparing the probe assembly

The ipsi probe tip and tubing are attached to the probe box at the factory. Connect the contra insert phone cable to the jack on the back of the base if contra reflex testing is to be performed. If contra testing will not be performed on a regular basis, it is not necessary to keep the contra phone attached to the system at all times.



- P1 - Green lamp:** *Blinking green* - Ready to start test.  
*Steady green* - Test in progress.
- P2 - Orange:** *Blinking orange* - Pressure leak.  
*Steady orange* - Occlusion.
- P1 and P2 off:** Test is finished.

### **WARNING**

Be sure that the Insert phone cable plug is inserted all the way into the jack to ensure proper operation.



**! WARNING**

To ensure the accuracy of calibration, the tygon tubing supplied with the ipsi probe assembly should not be cut or altered in any way. The system has been specifically calibrated to meet specifications with the tubing length supplied with your unit. A spare set of tubing is provided. If the replacement tubing supplied with your instrument is used, recalibration is unnecessary.

For the GSI 39 Combo probe, do not use the Flat probe eartips supplied with the GSI 39 handheld probe. The Contra phone assembly eartips and the probe eartips are the same and should be the rounded flange type.



**! WARNING**

A Cardinal Health/GSI provided Probe Tip must be used. Using the probe without the Probe Tip could result in injury to the subject.

## Front Panel Controls and Indicators



































Figure 2: Front panel

Legend / Label	Button	Description
<b>F1 / Print Screen</b>		button used to print the currently displayed page of memory or active test screen.
<b>F2 / Print All Memory</b>		Used to print all pages of data from memory.
<b>F3 / Paper Advance</b>		Causes paper to feed through printer; may be used to load paper or to provide space between printouts.
<b>F4 / PAGE</b>		Enters Page Mode: Pressing <b>F13</b> and <b>F14</b> scrolls through the test results stored in memory.
<b>F5 / M -</b>		Erases currently displayed page of data from memory.
<b>F6 / M - -</b>		Erases all pages of data from memory.
<b>F7 / Data Transfer</b>		Transfers test results to an attached computer.
<b>F8 / M+</b>		Save button; during <b>Audiometry</b> mode, causes the threshold information per frequency to be saved on the display; during <b>Program</b> mode, causes option to be selected.
<b>F9 / +10 dB</b>		Used to temporarily extend the intensity range by 10 dB; causes a large + sign to appear on the display indicating that the extended range has been selected.


Legend / Label	Button	Description
<b>F10 / Aud(iometry)</b>		Selects Audiometry mode.(Available in Versions 3 & 4 only). When in Audiometry, this button starts the Auto HL when held for 3 seconds.
<b>F11 / Headphone</b>		Selects the TDH39 calibration files for transducers. When the  button is pressed, the display will flash to ensure that you want to change the transducer selection. You must press the  button again to engage the TDH 39 headphone calibration file. The symbol  is shown on the center of the display if selected.
<b>F12 / Insert</b>		Selects the insert earphone calibration file for transducers. When the  button is pressed, the display will flash to ensure that you want to change the transducer selection. You must press the  button again to engage the insert earphone calibration file. The symbol  is shown on the center of the display if selected.
<b>F13 and F14 / Decrease and Increase Frequency</b>		Selecting  causes the cursor to move to the next lower frequency; selecting  causes the cursor to move to the next higher frequency.
<b>F15 / Steady</b>		Used during Audiometry mode to select a continuous test tone when Present Bar is depressed; causes the steady symbol to appear on the display.
<b>F16 / Pulsed</b>		Used during Audiometry mode to select a pulsed tone when the Present Bar is depressed; causes the pulsed symbol to appear on the display.
<b>F17 / FM</b>		Used during Audiometry mode to select a frequency modulated test tone when the Present Bar is depressed; causes the letters FM to appear on the display when selected.
<b>F18 / R</b>		Used to identify right ear is under test so that data stored in memory and/or printed is properly identified; for Versions 3 and 4, used to select right earphone for audiometry. An R will appear on the LCD.

<b>Legend / Label</b>	<b>Button</b>	<b>Description</b>
<b>F19 / L</b>		Used to identify left ear is under test so that data stored in memory and/or printed is properly identified; for Versions 3 and 4, used to select left earphone for audiometry. An L will appear on the LCD.
<b>F20 / Attenuator Knob (dB HL)</b>		Used to increase or decrease the intensity of the test tone presented in Audiometry mode; counterclockwise rotation causes the intensity to be lowered; clockwise rotation causes the intensity to be increased.
<b>F21 / Present Bar</b>		In Audiometry mode, press to present test signal to appropriate earphone; release to turn test tone off.
<b>F22 / TYMP</b>		Selects Tympanometry only mode.
<b>F23 / Tymp Reflex</b>		Selects Tympanometry and Reflex mode.
<b>F24 / 226Hz</b>		Selects 226 Hz for Probe Tone Frequency.
<b>F25 / 1KHz</b>		Selects 1000 Hz for Probe Tone Frequency.
<b>F26 / IPSI</b>		Selects an ipsilateral reflex test.
<b>F27 / CONTRA</b>		Selects a contralateral reflex test (available with Versions 2 and 3 only).
<b>F28 / Prog(ram)</b>		Selects Program mode screen which lists settings available for reflex presentation format, printout header format, audiogram vs. tabular format, display normal box, and identify frequency range for Audiometry mode.
<b>F29 / 500</b>		Selects 500 Hz as a stimulus during reflex testing.
<b>F30 / 1000</b>		Selects 1000 Hz as a stimulus during reflex testing. (Not available with 1000 Hz probe tone.)
<b>F31 / 2000</b>		Selects 2000 Hz as a stimulus during reflex testing.
<b>F32 / 4000</b>		Selects 4000 Hz as a stimulus during reflex testing.


Figures 3 through 7 show the individual display format for each test mode.



**Figure 3: 226 Hz Tympanometry screen.**




**Figure 4: 226 Hz Tympanometry/Reflex screen.**



**Figure 5: 1000 Hz Tympanometry screen.**

\* Margolis et. al.





**Figure 7: Audiometry screen.**



---

## Tympanometry testing information

As mentioned, it is a good idea to perform a test on a normal ear each day to make certain that your instrument is functioning properly. See **Biological Check** in Chapter 5 for details.

### Helpful hints

Tympanometry and acoustic reflex testing can be performed at any age. However, the technique used will vary with age. From three years through adult, tympanometry can be performed with little difficulty due to the cooperative nature of this age group. With the under-three-year population, a bit of ingenuity is required to keep the patient relatively quiet during the seconds required for the test. In all cases, distraction is the key to success. Anything that provides a sound and/or visual distraction should work.

Sucking on a pacifier or a bottle will help with the younger population. However, the tympanogram tracing will not appear as smooth due to the movement artifact. Having a parent hold an infant during testing will also help. For the 1000 Hz probe tone on infants, we recommend turning the **Auto Start** option off (factory default setting). This will allow you to position the probe and run repeated tests without removing the probe.

The key to success in all cases is to make sure that you are at eye level with the ear canal. Keep your hand steady and your eyes on the ear canal and probe lights until the test is over. It is a good idea upon first receiving your instrument to practice on a cooperative patient to gain confidence in its use. Once you feel comfortable with the probe, you are ready to handle any situation. Remain calm and success will follow.

## Obtaining a seal



**WARNING** A Cardinal Health/GSI provided Probe Tip must be used. Using the probe without the Probe Tip could result in injury to the subject.


Six different size eartips are provided with your instrument. The size eartip will vary with skeletal size of the individual it is to be used on. Generally speaking, the following criteria applies:

- Premie - 8 mm
- Newborn - 8 mm, 11 mm
- Pre-school - 11 mm, 13 mm
- School age - 11 mm, 13 mm, 15 mm
- Adult - 15 mm, 17 mm, 19 mm

**NOTE:** Before attempting to seal the entrance of the ear canal, visually inspect the opening to make sure that the canal is free of any obstruction. If the canal is completely plugged at the entrance or if fluid is running from the ear canal, tympanometry should not be attempted until the condition is cleared.


**NOTE:**Damage to the probe can result if fluid is pulled up into the probe with negative pressure.

1. Slip the appropriate size eartip onto the nose cone of the probe, making sure the rounded tip of the eartip sits flush with the tip of the nose cone (Figure 8).



**Figure 8: Positioning the eartip (226 Hz only probe)**

2. Move any hair away from the ear and pull upward and back on the pinna of the ear on an adult (pull downward and back on the pinna of a young child.) This tends to straighten out the ear canal and ensure better results. Keep the pinna in this position throughout the test sequence.
3. Make sure that the **green** lamp on the probe is blinking.
4. Position the probe up against the entrance of the ear canal, applying a gentle pressure to maintain a tight seal (Figure 9).




**Figure 9: Positioning the probe (226 Hz only probe).**

5. Watch the probe lamp. As soon as a good seal is obtained, the blinking green lamp will change to a steady glow and remain steady while the test is in progress.
6. Once the test sequence is over, all lamps on the probe will be turned off and the test result can be viewed on the instrument display before printing it out. It is now possible to remove the probe from the ear canal.
7. Note that the green lamp is now blinking again, signifying that you can run another test. Should you run into difficulty during the test, the probe lamps will inform you of a problem as follows:
  - **Green lamp:** Still blinking - seal has not been obtained to initiate the test sequence.
  - **Orange lamp:** The ear canal is not properly sealed and a large pressure leak exists.
  - **Yellow lamp:** The probe tip is occluded with cerumen or you are pressing the tip of the probe against the ear canal wall causing an occlusion.

In all cases, it is best to remove the probe, examine the tip for cerumen and clean it if necessary. A change of eartip size may also be appropriate. Start the test again.

## Combo Probe Insertion

Select the correct size eartip and position it fully on the probe. The eartip should be pushed firmly onto the tip of the probe until it is fully seated. The three probe tubes should be nearly flush with the top surface of the eartip.




**Figure 10: Inserting the Combo Probe.**

Insert the probe tip securely into the ear canal with a back-and-forth twisting motion. Pull the pinna upward and backward for adults, and down and back for children.

The probe tip should sit firmly within the ear canal without being held. If leakage occurs, a different size eartip may be needed.

## Audiometry testing information (Versions 3 and 4)

Prior to testing, ensure that the earphone cords are plugged into their appropriate connectors on the rear panel of the instrument. Both Headphones and Insert phones are available to use on the instrument. Select the appropriate transducer and the desired tone type (i.e., pulsed, steady, or FM).

 **CAUTION** Always handle earphones with care. Neither drop them nor permit them to be squeezed together. Severe mechanical shock may change their operating characteristics and require their replacement. Insert the earphone cords between the earphone cushions during storage to prevent damage from mechanical shock.


### Instructing the patient/subject

You should put the patient/subject as much at ease as possible before the test begins. In addition, it is important to try to make them understand how the test is to be conducted and what they will hear. For sake of uniformity, an unvarying explanation is advisable, for example:

“I am going to place these earphones over your ears. You will hear a variety of tones - some high, some low, some loud, and some very soft. Whenever you hear, or think you hear one of these sounds, raise your hand and then lower your hand when you no longer hear the sound.”

“Remember that although some of the tones will be easy to hear, others will be very faint. Therefore, you should listen very carefully and raise your hand even if the noise is very tiny.”

**NOTE:**Modify the instructions accordingly if Insert phones are being used or if indicating the sound is heard using the available handswitch.

 **WARNING** Any program aimed at obtaining reliable records of hearing thresholds should be staffed and supervised by appropriately-trained individuals. Training courses leading to certification are available for audiometric technicians in most urban areas.

### **Placement of earphones**

The most important thing to remember is that a good seal is required between the earphone cushion and the subject's/patient's head and ears. To increase the likelihood of a good seal:

- a. Eliminate all obstruction between the earphones and the ears (e.g., hair, eyeglasses, earrings, hearing aids, etc.).
- b. Adjust the headband so that it rests solidly on the crown of the subject's head and exerts firm pressure on both ears.
- c. Center the earphones carefully over both ears. The earphone with the red connector goes on the right ear. Take care to eliminate any visible gaps between the earphone cushions and portions of the individual's head and the ear on which the cushion rests.

### **Placement of Insert earphones**

1. Examine the ear canal for obstruction or excessive cerumen.
2. Make sure the sound tube is not blocked.
3. Place the black tubing of an ER-3A foam eartip completely onto the connector of the sound tube.
4. Roll the foam tip into the smallest diameter possible.
5. Insert the eartip well into the ear canal. Interaural attenuation is improved with deep insertion.
6. Allow foam to expand to acoustically seal ear canal.
7. Discard foam eartips after a single use.

**NOTE:** If using insert phones, ensure the appropriately sized foam tip is selected.

### **Response handswitch (optional accessory)**


If the optional handswitch is to be used, make sure that the handswitch connector is properly inserted into the jack on the rear panel. The instrument will display an appropriate symbol whenever the handswitch is pressed.

## Tympanometry/Reflex Test Sequence

This section describes the test sequences for all modes of operation. Since there are five versions of the GSI 39, all of the test sequences described may not apply to your particular unit. Section a. *Tympanometry Only Mode* applies to all five versions. Four of the five versions have tympanometry and ipsilateral reflex capabilities so the following three topics (*Tympanometry only mode*, *Tympanometry and Ipsilateral reflex*, and *Programming Ipsilateral acoustic reflex test frequencies* apply to Versions 1-4 units). The two topics later in this chapter (*Tympanometry and contralateral reflex* and *Tympanometry and ipsilateral/contralateral reflex*) are only applicable for individuals with Versions 2 and 3.

If a test sequence is not available on your system, “Invalid” will display on the LCD and those buttons will not activate the sequence. All systems are capable of being upgraded to add test modalities. Contact your local Cardinal Health/GSI Representative for more details on upgrade packages.

### a. Tympanometry only mode

1. Select the **Tympanometry only** mode by pressing  on the front panel. The display will immediately show the format for the tympanogram along with the summary information headers ECV, cm<sup>3</sup>, daPa, and GR. For 226 Hz Probe tone, the default scale for compliance is 1.5 cm<sup>3</sup>. If a compliance peak greater than 1.5 cm<sup>3</sup> is measured, the instrument automatically scales the compliance axis to 3.0 cm<sup>3</sup> so that more of the tympanogram data can be seen. For 1000 Hz Probe tone, the default scale is Baseline Off and will display a 5 mmho scale. If a larger compliance peak is measured, the instrument automatically scales the compliance axis to 10 mmhos.
2. Determine which ear is to be tested and select the appropriate ear (**R** or **L**) button so that the test results will be labeled properly. It is not possible to change the test ear once the probe is placed in the ear canal.
3. Examine the ear canal to determine the appropriate size eartip for the test and position the eartip on the probe. Be certain that the eartip is pushed as far down the probe tip as possible so that the eartip is flush with the tip of the probe. Position yourself so that you are at eye level with the test ear.
4. Note that the green lamp is blinking, which indicates that the instrument is ready to begin the test.
5. Place the probe up against the entrance of the ear canal so that it's opening is completely covered with the eartip and no visible leaks are apparent.



6. For 226 Hz, the test sequence begins once the instrument determines that a volume between 0.2 cm<sup>3</sup> and 5.0 cm<sup>3</sup> is present. This is indicated by the green lamp changing from blinking to a steady state. From this point on, hold the probe securely in this same position without any hand motion. Keep your eyes on the probe and the individual's ear.


At the start of the test, the pressure system establishes a pressure of +200 daPa within the ear canal. When this pressure is achieved, the instrument makes a measurement of ear canal volume. This information is valuable as it indicates whether a good seal has been established and helps to differentiate between two similar Tympanogram types (i.e., a fluid-filled middle-ear system and a perforated tympanic membrane). After the ear canal volume (ECV) is obtained, this compliance value is subtracted from the remaining compliance measurements so that a direct reading of the tympanogram compliance peak is possible.

The pressure sweep begins at the starting pressure of +200 daPa and proceeds in the negative direction at a rate of 600 daPa/second. Measurements of compliance are made continuously as the pressure sweep continues in the negative direction. The slope of the tympanogram increases as the measurement approaches the compliance peak. This signals the instrument to slow down the rate of pressure sweep to 200 daPa/second to ensure a more accurate reading of the compliance peak. After the peak compliance and pressure values are detected and stored, the tympanogram dips downward toward the baseline (i.e., 0 cm<sup>3</sup>) and the pressure sweep rate increases back to 600 daPa/second. The tymp sweep ends automatically when the compliance value returns to baseline and the pressure is at least -100 daPa. Only when the middle-ear pressure is very negative is it necessary for the pressure sweep to continue all the way down to -400 daPa. This automatic stop when the tymp compliance returns to baseline eliminates unnecessary pressurization of the ear and shortens the test time.

When the tympanogram is completed and the test is finished, the steady green lamp turns off and the tympanogram results are displayed.

For the 1000 Hz probe tone, the measurement taken at +200 daPa will be identified as C1 and will not be a calculated volume, but will be a mmho measurement.

For the 1000 Hz probe tone, the default setting for Auto Start is turned off. This setting is recommended to ensure the probe is placed properly in the ear canal prior to testing. To


begin the pressure sweep, press the  button.


The pressure sweep rate of 200 daPa/second remains steady through the peak measurement and will continue to -400 daPa. The green lamp will turn off when the test is completed.


7. It is now possible to remove the probe from the ear and to view the test results on the display.

The test results are stored automatically in a page of memory. The actual memory location number is determined by the number of tests that preceded this current test. For example, if this is the first test to be stored in memory, it will be assigned the number M1. If it is the third test to be stored in memory, it will be numbered M3, etc.

In addition to the tympanogram tracing, the screen displays the test summary information. For 226 Hz probe tones, this data includes the ear canal volume (ECV), the compliance peak in  $\text{cm}^3$ , the pressure at the peak of the tympanogram in daPa, and the gradient (GR) as a peak width value. This test result can be printed out immediately as a single test by

selecting the Print Screen Only  button or other tests can be run and saved before all

tests in memory are printed via the Print All  button. For 1000 Hz probe tones, the data accompanying the tympanogram will include a C1 value (in mmhos), a mmho compensated peak value (Peak - C1 = compensated peak in mmhos), and Peak Pressure in daPas.

**NOTE:** If a second tympanogram needs to be performed for 226 Hz, remove the probe and reinsert the probe. If Autostart is turned off on the 1000 Hz test, the probe does not need to be removed to run a second tympanogram. Press the  button to start a 2nd measurement. If the pump cannot run a 2nd tympanogram due to pressure equalization, **remove** will appear on the bottom right of the LCD. Remove the probe and reinsert for the next tympanogram.

## b. Tympanometry and Ipsilateral Reflex

The default parameters for this test are tympanometry followed by an ipsilateral acoustic reflex test at 1000 Hz (2000 Hz for 1000 Hz probe tone).

Once a seal is obtained, the tympanometry sequence is initiated. (See topic *Tympanometry only mode* earlier in this chapter for details). As long as no large leak is encountered during tympanometry (orange lamp illuminated) and no occlusion is detected (yellow lamp illuminated), the test automatically sequences on to the reflex portion of the test as follows.

- a. For 226 Hz probe tones, the pressure from the tympanogram peak compliance is re-established within the ear canal and is offset by -20 daPa so as to avoid any problems with extremely sharp tympanogram slopes (+20 daPa for positive peak pressure). For 1000 Hz probe tones, the pressure is re-established to 0 daPa for reflex measurements.
- b. With the air pressure held constant throughout the reflex test sequence, the lowest intensity level for the starting frequency is presented and a measurement of compliance change is made. If the compliance change of at least 0.05 cm<sup>3</sup> for the 226 Hz probe tone and 0.1 mmho for the 1000 Hz probe tone is measured, this reflex intensity level is stored in memory.
- c. If no other frequencies were selected for the test, the Tymp Reflex sequence ends here. The green lamp is no longer illuminated indicating that it is time to remove the probe from the ear. The display will indicate the reflex test result as a Yes, as an HL value, or as an HL value plus a tracing of the reflex response curve. The default setting established in the Program mode determines the manner in which the reflex result is displayed. See *Program Mode* in this chapter.

- d. If no response is measured (i.e., for 226 Hz probe tone, a compliance change of at least 0.05 cm<sup>3</sup> was not detected) at this lowest intensity level, the intensity level of the stimulus is automatically increased by 10 dB. During this second presentation, a measurement of compliance change is also made. If a response is detected, the test sequence for this frequency ends and either the result is displayed on the screen or the test proceeds on to the next frequency selected. However, if once again no response is detectable, the intensity level is increased by 10 dB (e.g., 1000 Hz Ipsi = 105 dB HL) and the stimulus is presented.
- e. After the compliance measurement is made and a response is detected, the intensity level is stored as the reflex test result and displayed on the screen. If no response is detectable at the highest intensity level, either a No or an NR (depending upon the Program mode setting) is indicated on the screen next to the frequency tested label. If during any of the three stimulus presentations a large pressure leak develops, an NT will appear on the screen next to the reflex test frequency and the test sequence is aborted.
- f. The same sequence is followed for each test stimulus selected.

**NOTE:** To change the test default frequencies, see *226 Hz Reflex* and *1000 Hz Reflex* in *Program Mode*, Chapter 4 of this manual.

The intensity levels available vary with the frequency selected ipsilaterally as follows:

IPSI	Intensity Levels
<b>500 Hz</b>	80, 90, 100 dB HL
<b>1000 Hz</b>	85, 95, 105 dB HL
<b>2000 Hz</b>	85, 95, 105 dB HL
<b>* 4000 Hz</b>	80, 90, 100 dB HL

**226 Hz Probe Tone**

IPSI	Intensity Levels
<b>500 Hz</b>	80, 90 dB HL
<b>2000 Hz</b>	85, 95 dB HL
<b>4000 Hz</b>	80, 90 dB HL

**1000 Hz Probe Tone**




**NOTE:**1000 Hz is not available when using the 1000 Hz Probe Tone option.

**NOTE:**Although four frequencies are available during the tympanometry and ipsilateral reflex test mode, most situations require only one or two frequencies to be tested. You can choose from a selection of the most commonly used frequencies. However, it is strongly recommended that you select only one to two frequencies per test. Holding the probe in the same position for the length of time it takes to test four frequencies may become a problem for both you and the individual being tested.

\* Only 80 and 90 dB HL are available on the Combo Probe.

## c. Temporary programming of ipsilateral acoustic reflex test frequencies

As mentioned earlier, the instrument defaults to a 1000 Hz Ipsilateral test stimulus when the **TYMP REFLEX** button is first pressed after receipt from the factory. However, any combination of the four available frequencies (500, 1000, 2000, 4000 Hz) for 226 Hz Probe Tone and three frequencies (500, 2000, 4000 Hz) for 1000 Hz Probe Tone can be selected either temporarily or as revised default parameters. To temporarily modify the default condition:

1. Press the **Tymp Reflex**  button.
2. Select the test **frequencies** by pressing the desired Frequency button (e.g.,  or ). Pressing the Frequency button a second time will deselect that frequency from the test sequence. Test frequencies must be selected before the probe is in the ear.






Each frequency selected will be indicated on the display. For example, if 2000 Hz is selected along with 1000 Hz, the label “I 1000” will appear at the top of the first column of numbers for reflex and “I 2000” will appear directly below it. If 500 is also selected, the screen will be modified so that “I 500” appears at the top of the first column of reflex numbers, “I 1000” will appear directly below “I 500” and “I 2000” will appear at the top of the second column of reflex numbers and directly to the right of “I 500” and so on.

**NOTE:**There is a pattern to the way in which these frequencies are positioned on the screen; the lowest frequency will be placed at the top of the first column for reflex results followed by the next lowest frequency. If more than two frequencies are selected, the third and fourth frequencies will be placed in the second column for reflex results in a low to higher frequency order.

To change the default setting, see the *Program Mode* in Chapter 4 of this manual.

#### d. Tympanometry and Contralateral Reflex (Version 2 and 3)

To select tympanometry and contralateral reflex testing:

1. Press the **Tymp Reflex**  mode button. This initializes the GSI 39 to perform a tympanogram along with reflex measurements. The factory default setting for reflexes is 1000 Hz Ipsilateral presentation.
2. To temporarily change the system to perform Contralateral Reflexes only, you must first deselect the 1000 Hz Ipsi. To deselect the 1000 Hz IPSI, simply press the  button, "I 1000" should no longer be displayed on the right side of the LCD.
3. Press the **CONTRA**  button. This causes the letter **C** to appear in front of the frequency labels.
4. Select the test frequencies by pressing the desired Frequency button (e.g.,  or ). Pressing the Frequency button a second time will deselect that frequency from the test sequence. Test frequencies must be selected before the probe is in the ear canal.
5. Before initiating this test sequence, select the appropriate **size eartip** for the contralateral insert phone from the color-coded eartip container. The size selected should be such that the insert phone can be tightly fitted into the ear canal.
6. Push the selected eartip firmly onto the insert phone. Be sure to carefully position the insert phone within the ear canal as the calibration depends upon a proper seal of the ear canal.
7. Select the test ear by pressing **R** or **L**. According to general convention for recording contralateral reflexes, the test ear is the ear where the probe is positioned and the stimulus ear is the ear that contains the contralateral insert phone. If the contralateral insert phone is placed in the left ear, the test ear is the right ear since this is the ear from which the reflex response is to be measured. Press the button that corresponds to the ear with the probe to select the test ear.

To initiate the test:

1. Position the insert phone securely within the ear canal to receive the contralateral reflex stimulus.
2. Position the probe up against the ear canal of the ear under test. Note that the green lamp changes from a blinking to a steady state once the test begins.
3. Keep your eye on the probe and the ear canal throughout the test sequence. The test begins with the tympanogram and is followed immediately by the contralateral acoustic reflex test.

For **226 Hz Probe Tone**, the pressure value used within the test ear throughout the contralateral stimulus presentations is the peak pressure obtained during the tympanogram offset by -20 daPa (+20 daPa if peak pressure is positive).

For **226 Hz Probe Tone**, as with ipsilateral reflex testing, a compliance change of 0.05 cm<sup>3</sup> indicates a reflex response is detected. Up to three intensity levels per frequency selected are presented. Recall that it is only necessary to go beyond the first and lowest level if a change of 0.05 cm<sup>3</sup> cannot be detected during the stimulus presentation period. The format in which the test results are displayed on the screen is determined by the default setting chosen in the Program mode (i.e., yes/no, dB HL, or dB HL and response curve). The three intensity levels available per frequency are the same for all four (500, 1000, 2000, and 4000 Hz) possible frequencies i.e., 90, 100, and 110 dB HL.

For **1000 Hz Probe Tone**, the reflexes will be measured at 0 daPa and a compliance change of  $\geq 0.1$  mmho is required. Up to 2 intensity levels per frequency are presented.

**NOTE:**Remember that the second or third intensity level presentations occur only if a response is not detected at the prior intensity level. The test is over once the green lamp on the probe is no longer illuminated.



**e. Tympanometry and ipsilateral/contralateral reflexes (versions 2 and 3 only)**

As described earlier, this test sequence can be selected either temporarily or set as the default sequence. If both ipsilateral and contralateral testing are performed only with certain patients, it is advisable to only change the test parameters temporarily on an as needed basis. However, if your test protocol calls for ipsilateral and contralateral testing with all patients, it is advisable to change the default settings. (See *226 Hz Reflex* or *1000 Reflex* under **Program Mode** in Chapter 4).


**Ipsi and Contra Acoustic Reflex testing**



There are 4 frequencies of either Ipsi or Contra stimulus presentations available. They can be all Ipsi, all Contra or a combination of both Ipsi and Contra. Ipsi will always be presented first and the frequencies will always go from low to high. Any combination of Ipsi and Contra frequencies can be programmed into the 4 stimulus conditions. For example, you could have the following:

I	500	C	1000
C	500	C	2000

or

I	500	I	2000
I	1000	C	1000

Press the  button to select the Tymp/reflex mode.

Press the  button and select and deselect ipsi reflexes using the frequency buttons. Press the  button and select contra reflexes required using the frequency buttons.

**NOTE:** Only 4 reflexes are allowed. It is not possible to select a 5th reflex. To choose different reflexes, you must first deselect those reflexes not desired.

Before initiating the test:

Position the insert phone securely within the ear canal to receive the contralateral reflex stimulus.

Position the probe up against the ear canal to receive the ipsilateral stimulus.

Once the green lamp changes from a blinking to a steady state, the test sequence begins. First, a tympanogram is obtained and then, for 226 Hz pulse tone, the peak pressure from the compliance peak offsets by -20 daPa (or a + 20daPa for a positive pressure peak) is re-established within the ipsilateral ear canal (i.e., ear containing the probe). For 1000 Hz, the system is reset to 0 daPa. The reflex sequence begins automatically by starting with the lowest ipsilateral test frequency and is followed by a second ipsilateral test frequency if selected. After the ipsilateral reflex tests are completed, the instrument sequences automatically to the contralateral reflex test stimuli. Once again, the lowest frequency is presented first and is then followed by the next frequency. Keep your eye on the ear canal where the probe is positioned. Once the green probe lamp is no longer illuminated, the test is complete and it is possible to remove the probe and the insert phone from their respective ear canals. The reflex test results can now be observed on the display screen. The format in which the ipsi and contra reflex test results are displayed is dependent upon the setting chosen in the Program mode.

### **Exit tympanometry/reflex**

To exit **Tymp Only Mode**:

Select **Tymp Reflex** or **Audiometry Mode**. Note that the appropriate screen appears on the display.

To exit **Tymp/Reflex Mode**:

Select **Tymp** or **Audiometry Mode**. Note that the appropriate screen appears on the display.





---

## Audiometry Test Sequence (Versions 3 and 4 only)

### To enter the Audiometry mode

1. Press the **AUD** button. Note that the display changes from a Tymp or Tymp/Reflex format to an audiogram format.

### Transducer Selection

Select the transducer to be used for the Audiometric testing. Press  to select Headphones or  to select Insert Phones. The LCD will flash a picture of the transducer choice until the transducer button is pressed a second time. With one set of output jacks for the transducers, two buttons allow separate calibration files to be accessed. Be sure that the transducers that are connected to the back of the GSI 39 are the same as the selected transducer from the front panel. If Headphones are selected, a  will appear in the center of the LCD. If Insert Phones are selected, a  will appear in the center of the LCD.

The settings for the frequencies available during audiometry are set in the Program mode as 125 through 8000 Hz (normal) or 500 through 6000 Hz (narrow). The factory default setting is the normal frequency range of 125 through 8000 Hz. Upon entering the audiometry mode, the starting frequency is automatically selected to be a steady signal of 1000 Hz at 0 dB HL.

## **GSI 39**

---


You can change the signal format temporarily from steady (continuous) to a pulsed or frequency modulated tone. These alternative tone formats remain selected as long as you remain within that audiometric test. Once you leave that specific test by selecting either Tymp or Tymp Reflex and re-enter the audiometry mode, the tone type returns to steady. The display indicates a **continuous bar** when steady is selected, a **dashed bar** when pulsed is selected, and the letters **FM** when frequency modulation is selected.



The audiometry defaults to testing the right ear first. To start with the left ear, it is necessary to press the **L** button after entering the audiometry mode. Since the audiometry mode defaults to 1000 Hz at 0 dB HL, the cursor is positioned at the corresponding location on the audiogram.



Please note that even though you may have selected the tabular format for the audiometric test results on the printout, the screen always appears in the audiogram format.

## To change the frequency

1. Press the  Hz button.

If the  Hz button is pressed once momentarily, the frequency increases to the next frequency in the range.

If the  Hz button is held down continuously, it is possible to quickly scroll through the available frequencies. Note that if the button is held down past the 8000 Hz in the normal range (6000 Hz for the narrow range), the frequency scroll wraps around to the lowest frequencies (i.e., 125 Hz with the normal range and 500 Hz with the narrow frequency range). The reverse occurs if the  Hz button is pressed.

In addition to changing the frequency, the  and  buttons change the position of the cursor on the audiogram. The frequency value of the cursor position on the audiogram is displayed on the right side of the screen.



## To change the intensity level of the test tone





1. Rotate the **dB HL** knob in the clockwise direction to increase the intensity level in 5 dB steps; rotate the knob in the counterclockwise direction to decrease the intensity level in 5 dB steps.


The cursor on the audiogram moves up and down accordingly. Also, the dB level displayed above the frequency value on the right-hand side of the audiogram changes.

For each frequency, there is a fixed intensity range available while rotating the **dB HL** knob as follows:

Frequency	Intensity Range
125 Hz	-10 to 50 dB HL
250 Hz	-10 to 70 dB HL
500 to 6000 Hz	-10 to 90 dB HL
8000 Hz	-10 to 70 dB HL




It is possible to extend the intensity range per frequency by 10 dB by pressing the  button. The button may only be selected when the intensity level is set to the highest value in the normal range. For example, with the test tone of 1000 Hz, the normal intensity limit is 90 dB HL. When the intensity knob is rotated clockwise to select beyond 90 dB HL, the intensity value above the 1000 Hz to the right of the audiogram flashes indicating the maximum intensity limit has been reached. To go beyond 90 dB HL, select the  button. A large **+** sign appears on the screen below the 1000 Hz value. The **dB HL** knob can be rotated through two additional positions, 95 and 100 dB HL. Rotating the **dB HL** knob to the next position beyond 100 dB causes the intensity value 100 flashes on the screen to the right of the audiogram; this indicates that the maximum dB HL for the extended range has been reached. If the dB HL is rotated one more position beyond the flashing 100 dB position, the letters NR appear next to the letters dB above the 1000 Hz. This permits the selection of the no response (NR) symbol on the audiogram during testing. The extended range remains selected until either the intensity level for that particular frequency (e.g., 1000 Hz) is brought down 5 positions below the maximum dB HL value (e.g., 65 dB HL for 1000 Hz) or the frequency is changed.

To save the threshold for a frequency, press the  button. The appropriate symbol (**0** for right ear and **X** for left ear) for the test ear will replace the cursor. If no response (NR) was detectable over the intensity range available, an arrow is attached to the 0 or X symbol on the audiogram. It is possible to repeat a threshold check for any frequency by returning to that frequency by way of the  and  **Hz** buttons. The last threshold obtained and saved with the  button becomes the value saved in memory and is the value printed out on the audiometric test results.

To present the tone to the test ear, press the **Present** bar. A speaker symbol  appears in the center of the screen for as long as the **Present** bar is depressed.

**NOTE:** Although the printout will combine the right and left ear test results on the same audiogram or table, the screen can display only the results from one ear at a time. Therefore, if an ear button (**R** or **L**) is selected while you are still testing a particular ear, the screen will change to a new audiogram. Once this happens, it is not possible to return to an incomplete audiogram to complete the test sequence.

## Screening audiometry

1. Carefully position the earphones over the individual's ears so that the **red phone** covers the right ear and the **blue phone** covers the left ear.
2. Be sure that nothing is obstructing each earphone such as earrings, eye glasses or a hearing aid.
3. Instruct the person being tested to raise a hand or a finger (or press the optional **Handswitch**) whenever a tone is heard.
4. Encourage the patient to respond even if he/she is not sure whether a tone is heard.
5. Select the ear to be tested with the **R** (right) or **L** (left) button.
6. Select the desired screening intensity by rotating the **dB HL** knob to the appropriate position. The American Speech Language and Hearing Association recommends 20 dB as the screening level for school-age children.
7. Select the starting frequency by pressing the  or  **Hz** buttons.
8. Present the tone by pressing the **Present** bar.
9. If the individual fails to respond, increase the intensity by 10 dB and try again. Press the  button at the intensity level where the individual responded.
10. Continue the procedure for all the desired frequencies.




## Audiometric Threshold

The GSI 39 offers two ways to collect Audiometric Threshold. The system can be used in a Manual mode or an Automatic Hearing Level mode (**Auto HL mode**). In the Manual mode, the intensity, frequency and presentation of the stimulus are controlled by the tester. In the Auto HL mode, the system presents stimuli based on responses from the Patient Response switch.

### Manual Threshold Audiometry

1. Carefully position the earphones and select the ear to be tested.
2. Familiarize the individual with the test procedure by presenting a tone of 40 dB HL at 1000 Hz.
3. Decrease the intensity in 10 dB steps until the person no longer responds or until you reach 0 dB HL.
4. When you believe the individual understands the procedure (i.e., raise your hand/finger when you hear a tone) proceed with the test procedure. The level of the first presentation generally is 10 dB below the level at which the individual responded during the familiarization procedure.
5. Starting at the desired test frequency, present the tone for a period of one or two seconds.

6. If a response is indicated,
  - a. Decrease the intensity of this same test frequency by 10 dB and present the tone again for one to two seconds.
  - b. If no response is indicated, increase the intensity by 5 dB. Present the tone again.
  - c. If no response is indicated, increase the intensity by another 5 dB.
  - d. If a response is indicated, this is the second time that the individual responded to the same intensity level. Repeat the sequence of down 10 dB and up in 5 dB increments to determine if a correct response is again detected at the same intensity level. The threshold is considered to be the minimum level at which a response has occurred two out of three times. Press the  button when this intensity level is indicated on the screen above the test frequency to signify that the threshold level for that frequency has been reached. Note that the appropriate symbol (**0** = right, **X** = left) appears at the correct intensity level where the threshold was determined.
7. Repeat this test sequence for each frequency to be tested.
8. Once the thresholds have been obtained for all the desired frequencies, select the other ear and repeat the sequence. Note that the display changes to a new screen for storing the second ear's results. The test protocol follows a down 10 dB and up 5 dB sequence to arrive at the threshold level.

## Automatic Hearing Level

The Automated Hearing Level Procedure (Auto HL) was designed to control the stimulus presentation via the patient control. The software determines the presentation level of the stimulus based on the Hughson-Westlake Threshold estimation procedure (reference). The patient is given a button that is hard wired to the GSI 39 and should be instructed to hold the button down when they hear the tone and release the button when the tone goes off. In this procedure the stimulus level is decreased 10 dB each time the patient presses the button and increased by 5 dB when the button is not pressed. The GSI 39 will present the stimulus and increase or decrease the intensity of the stimulus based on the patient response. The GSI 39 keeps track of the response/no response stimuli and determines the hearing threshold based on the data.

## Theory of Operation


The following bullet points describe the stimulus presentation patterns and patient response validity:

- The stimulus (tone) on time is fixed at 1.5 seconds.
- The Inter-stimulus interval is randomized between 3 and 5 seconds.
- When a valid response occurs, the intensity for the next stimulus presentation is decreased 10 dB from the current one. When no valid response occurs, the intensity of the next stimulus presentation is increased 5 dB. This is based on the Hughson-Westlake down 10, up 5 dB rule used by most audiologists during threshold testing.
- The system will determine the response to be valid if the patient response switch is pressed during the stimulus or for 2 seconds following the stimulus off time.
- The system will determine the patient response to be invalid based on the following occurrences:
  - a. The patient response switch is pressed during the stimulus on time, but not released before the start time for the next intensity presentation.
  - b. The patient response switch is pressed and released only during randomized inter-stimulus interval.
  - c. The patient response switch is pressed and released more than 2 times during the stimulus on and completion of inter-stimulus interval.

Threshold results are displayed as they are obtained for each frequency. Once the first ear test sequence is completed, the audiometric thresholds for all the frequencies tested are stored in memory. At the start of the second test ear sequence, the results on the LCD will be cleared to display the second ear results. Once the second ear sequence is completed, the entire audiogram containing thresholds for both ears are stored in memory.



The threshold series for any frequency will be considered invalid if a threshold is not achieved within 18 stimulus presentations, or if the retest result at 1000 Hz doesn't agree within 5 dB of the first result. If the threshold results are considered invalid the system will exit the Auto HL procedure. The audiogram results obtained thus far will be kept and displayed so that the test can be completed manually.

## Performing the Auto HL Procedure

1. Instruct the patient to press the button on the handswitch when they hear the tone and release the button when the tone goes off.
2. Carefully place the headphones or insert earphones.
3. To begin the Auto HL procedure, press the **AUD** button and hold it for 3 seconds. The words **Auto HL** will be displayed at the bottom right corner of the LCD indicating that the Auto HL procedure has been engaged. The first stimulus will be presented when the **AUD** button is released. When a signal is presented, the Speaker  icon will display on the LCD.
4. When all frequencies have been successfully tested the **Auto HL** will disappear from the LCD indicating the test is finished.

## Exit audiometry


There are two ways to exit the audiometry mode:






- a. Select the **Tymp** mode  button
- or -
- b. Select the **Tymp Reflex** mode .

For details on programming the Auto HL procedure, refer to *Programming the Auto HL Procedure* in Chapter 4.





---



## Tests in memory


The Tympanometry and Tympanometry Reflex test results are automatically stored in memory when the test sequence ends. Audiometric test results are stored in memory when  is pressed. A total of 12 memory pages are available with the GSI 39. Each Tympanometry, Tympanometry/Reflex or individual ear in audiometry is assigned a page in memory. They are labelled M1 - M12.



To review the individual test results, press the  button and enter "Page mode". The memory number is located in the upper right-hand corner of each screen. If, for example, only five tests were stored in memory, only five memory locations can be scanned. The memory can be scanned a page at a time by pressing the  or  button once and observing the result. The entire memory can be scrolled through by holding the  or  buttons down continuously.

### Memory erase

If there is a particular test result that you wish to delete before printing, enter the Page mode by pressing . Press  or  to display the test result and press . This erases that particular test result from memory. The LCD displays a blank screen for erased memories with the memory location number located at the top right corner. Upon exiting from the Page mode, the stored memories reshuffle and replace the empty memory with the remaining tests in the order in which they were obtained. The Page mode will be exited once you press the **PRINT**





**ALL**  or **ERASE ALL**  buttons or any button that would normally begin the setup of a new test. Page Mode is ready only. No changes can be made to audiometric results.





If you should wish to erase all tests from memory, press the **ERASE ALL**  button.

**NOTE:** Be certain that you wish to remove all tests from memory before pressing the  button because the erasure occurs immediately upon pressing the  button!

## Printing test results

The printout will begin with a header, if it is selected, in the program mode (i.e., GSI 39 or a custom header designed by you). The next two lines contain space for recording the individuals name and the test date. This is followed by the test results in the order that they were obtained/selected.

Either a single test can be printed from memory or the entire group of tests in memory can be printed. To print a single test from memory, use the **PAGE**  button to enter the Page mode and the  or  button to arrive at the desired test result to print. Once this test is displayed, press the **PRINT SCREEN**  button.

To print all tests in memory, press the **PRINT ALL**  button. When **PRINT ALL**  is pressed and two audiogram tests are stored in memory, they will combine under the following conditions. There must be one left test and one right test sequentially stored in memory. A left and right audiometric pair of tests will not be combined if they are separated in the memory by a Tympanometry test. Therefore, when tests are erased, the result could cause a change in (left, right) or (right, left) sequence with Audiometric tests. This would result in the wrong audiometric tests being combined when **PRINT ALL**  is selected. Prior to selecting **PRINT ALL**  you should scroll through the tests in memory to determine where the audiometric tests are located. This will help you avoid combining tests from different patients.

To avoid this issue, you may want to press **ERASE ALL**  before starting a new test patient.

# Chapter 4




# Program Mode








---

Blank page.



## Program Mode

To enter the program mode, press the **PROGRAM** button located on the front panel. There are two screens for the Program mode. To move to the second page, press the Increase Frequency  button or turn the Attenuator knob  until the cursor is next to the arrow on the bottom right column. Press  to enter Page 2.

Basic button functions for moving through the Program menu	
 (Cursor)   (Attenuator Knob)	Moves the cursor sequentially through the list of options on the screen.
	Toggles the option on or off. An asterisk (*) appears to the left of the item to denote the item has been selected. Pressing  again removes the asterisk, which deselects the item.
 (Page)	Use this button to move to the submenu or next page of a menu.
Save	The word Save should appear on the lower right corner of the LCD after the  button has been selected.
...	Indicates there is a submenu. Select  to enter the submenu.

## Program Mode Menu Items

The following screen appears the first time you enter the program mode.

PROGRAM MENU PAGE 1		PROGRAM MENU PAGE 2	
PROBE HZ...	AUD RANGE NORMAL	DATA XFER CONFIG...	INTERNAL PRINTER
TYMP OPTIONS...	AUD RANGE NARROW	POWER UP SETTINGS	EXTERNAL PRINTER
REFLEX DISPLAY...	PRINT - AUDIOGRAM	PRN HEADER GSL...	RESET TO DEFAULTS
226 HZ REFLEX...	PRINT - AUD TABLE	PRN HEADER OFF	→
1 kHz REFLEX...	DEF XDUCER TDH 39	PRN HEADER CUSTOM	
AUTO HL SETUP	DEF XDUCER INSERT		
LANGUAGE	→		


**NOTE:**Pushing the **Print** button while in Program mode will print out the currently programmed settings.

**NOTE:**The GSI 39 is sold in 5 versions, each containing different testing modalities. When navigating through the Program menu, features that are not available on the GSI 39 version you purchased will be represented by **invalid** on the LCD.

**NOTE:**Factory Default settings are listed at the end of this chapter.

## Program Menu Page 1 Option Descriptions

### PROBE HZ . . .

This submenu allows you to select the probe frequency that you want the GSI 39 to display when the system is first turned on. Move the cursor to the desired frequency and press the  button to save the selection.

**226 Hz**                      →

**1000 Hz**

### TYMP OPTIONS . . .


This submenu allows you to select tympanogram display and test options.

NORMAL BOX ASHA		BASELINE ON	1k
NORMAL BOX OFF		BASELINE OFF	1k
NEWBORN NRM ON	1k	AUTOSTART ON	1k
NEWBORN NRM OFF	1k	AUTOSTART OFF	1k
50th PERCNT ON	1k	→	
50th PERCNT OFF	1k		


### NORMAL BOX ASHA/NORMAL BOX OFF

For 226 Hz Probe Tone, it is possible to have the Normal Box, as defined by ASHA, appear on the tympanogram screen and printout. The boundaries for this Normal Box are -150 daPa to +100 daPa and 0.2 cm to 1.4 cm<sup>3</sup>.


**NOTE:** A compliance value of 1.5 cm<sup>3</sup> or greater will turn off the ASHA normal box automatically.

**NORMAL BOX ASHA** is the factory default setting. To select the **NORMAL BOX OFF**, move the cursor next to the selection and press the  button to save. **Saved** should appear on the bottom right corner of the LCD and an “\*” should appear next to the **NORMAL BOX OFF** option to denote the selection.

## **NEWBORN NRM ON 1k / NEWBORN NRM OFF 1k**

For 1000 Hz Probe Tone, it is possible to have the normal box, as described by Margolis, et.al., appear on the tympanogram screen and printout. The **NEWBORN NRM ON 1k** is the factory default. The **NEWBORN NRM ON 1k** are represented on the display by dashed lines at the 5th or 95th percentile. To select the **NEWBORN NRM OFF 1k** option, move the cursor next to the selection and press the  button to save. **Saved** should appear on the bottom right corner of the LCD and an “\*” should appear next to the **NEWBORN NRM OFF 1k** option to denote the selection.

## **50th PERCENT ON 1k / 50th PERCENT OFF 1k**

In the normative data put forth by Margolis, et. al. a dashed line representing the 50<sup>th</sup> percentile of the infant population is presented. This dashed line can be turned off by selecting **50<sup>th</sup> PERCENT OFF** and pressing the  button to save the selection.


## **BASELINE ON 1k / BASELINE OFF 1k**


**BASELINE ON 1k:** The C1 value in mmhos is obtained at +200 daPa and then subtracted from the tympanogram tracing so that it begins at the 0 mmhos position on the tymp display (compensated tymp data is displayed).

**NOTE:** This display invalidates the Newborn Norm display. When selecting **BASELINE ON 1k**, you must also select **NEWBORN NORM OFF 1k**.

**BASELINE OFF 1k:** The C1 value in mmhos is obtained at +200 daPa; the tympanogram tracing begins at this amplitude at the +200 daPa position (uncompensated tymp data is displayed).

**AUTOSTART ON 1k / AUTOSTART OFF 1k**

For the 1000 Hz option, you can turn off the Autostart option by selecting the **AUTOSTART OFF 1K** option. Move the cursor to **AUTOSTART OFF 1K** and press the  button to save the selection.

If **AUTOSTART** is **Off**, press the  button to start the test.

**NOTE:** Turning off the Autostart feature will allow you to situate the probe in the ear before the testing can begin. This feature also allows for repeat tympanograms to be obtained quickly without having to remove the probe.


**REFLEX DISPLAY**

Reflex test results can be displayed and printed in three different formats:

• **Reflex dB HL plus curve**

The default setting for this grouping is **Reflex dB HL plus curve**. All reflex test results will appear on the display and printout with the following information:


- a. **I** (Ipsi) or **C** (Contra) if available and selected
- b. **Frequency**: 500, 1000, 2000, or 4000 Hz
- c. **Intensity** level where response was detected
- d. Tracing of actual response curve.



**Figure 11: Display format for TYMP/REFLEX Test (Reflex test results with dB HL value and tracing).**

• **Reflex dB HL only**


If **Reflex dB HL only** is selected, the stimulus frequency, stimulus routing, and the dB HL level for the reflex will appear on the display and printout.



**Figure 12: Display format for TYMP/REFLEX Test (Reflex test results given in dB HL).**

• Reflex yes/no

If **Reflex yes/no** is selected, the dB HL result will be replaced with the word **yes** (response detected at one of three levels) or **no** (no response detected).



**Figure 13: Display format for TYMP/REFLEX Test (Reflex test results given as Yes or No).**

When the reflex test cannot be performed, due to a leak or early extraction of the probe, a “NT” will appear next to the frequency.

To select a different setting for reflex format:


1. While in the Program mode, move the cursor down to the setting that you wish to select for your own default criterion.
2. While the square cursor is positioned in front of your desired setting, press the  button.

The word **SAVED** appears in the lower right corner of the screen. The previous setting is deselected. An asterisk (\*) is displayed beside the new default setting.

## 226 Hz REFLEX

<b>Ipsi</b>	500	<b>Contra</b>	500
<b>Ipsi</b>	1000	<b>Contra</b>	1000
<b>Ipsi</b>	2000	<b>Contra</b>	2000
<b>Ipsi</b>	4000	<b>Contra</b>	4000
		→	

This feature allows you to select the stimuli and signal routing of the acoustic reflexes as default settings. To select the frequencies, move the cursor next to the selection and press the


 button to save. **Saved** will appear on the bottom right corner of the LCD. An \* will appear next to the selected stimulus routing and frequency. The system will allow 4 stimulus choices of any combination (i.e., ipsi or contra) for display and printing.

## 1000 Hz REFLEX

<b>Ipsi</b>	500	<b>Contra</b>	500
<b>Ipsi</b>	2000	<b>Contra</b>	2000
<b>Ipsi</b>	4000	<b>Contra</b>	4000
		→	

**NOTE:**1000 Hz reflex stimulus is not available for the 1000 Hz probe option.


This feature allows you to select the stimuli and signal routing of the acoustic reflexes as default settings. To select the frequencies, move the cursor next to the selection and press the

 button to save. **Saved** will appear on the bottom right corner of the LCD. An \* will appear next to the selected stimulus routing and frequency. The system will allow 4 stimulus choices of any combination (i.e., ipsi or contra) for display and printing.






## AUTO HL SETUP


### Programming the Auto HL Procedure

Navigate the cursor to the **Auto HL Setup** line located on the **Program Mode Screen 1** and press the  button. The following submenu will appear:




**Test Frequencies (Hz) . . .** →  
**Intensity Range (dB HL) . . .**  
**Start test ear . . .**  
**Scoring rule . . .**  
**Tone Format . . .**

Place the cursor next to the line item, press the  button to enter the submenu item. Auto HL features are selected in the Program mode by placing the cursor next to the parameter and pressing the  button to engage the selection.

At any time you wish to exit this submenu, move the cursor to → and press .

**Test Frequencies (Hz):** This submenu allows you to select the frequencies to be tested during the Auto HL procedure. Move the cursor to the frequency and press the  button to select or deselect the frequencies for presentation during the Auto HL procedure. An asterisk next to the frequency denotes that it has been selected for presentation. The submenu will appear as follows with the factory default settings:

Test frequencies (Hz)	
125	*2000
250	*3000
*500	*4000
750	*6000
*1000	8000
1500	Return to Auto HL Set




**Intensity Range (dB HL):** This submenu allows you to determine the minimum and maximum decibel level (HL) that will be presented during testing. To change the Min. dB (lowest level), place the cursor on that line and turn the **HL knob** on the front panel to the desired level. Press the  to move the cursor to the Max dB line and use the **HL knob** again to change the maximum level. Press the  button to move the cursor to the **Return to Auto HL Setup** and press  to exit the submenu. Asterisks on this menu denote factory default settings.

Intensity Range (dB HL)	
Min. dB: 0*	Return to Auto HL Set up
Max dB: 90*	

**NOTE:** Setting the **Min. dB range to 20** and the **Max dB range to 45** will allow you to perform a quick screening procedure using the Auto HL feature.



## LANGUAGE


Six language selections are available. Use  to highlight **LANGUAGE** and press the  button to enter the Language submenu and then move the cursor to the language you wish to select. Press the  button to activate the selected language. The following languages are available.

ENGLISH  
SPANISH  
ITALIAN  
GERMAN  
FRENCH  
PORTUGUESE



The new language will activate immediately upon exiting the Language submenu.

## AUD RANGE NORMAL/AUD RANGE NARROW

All eleven frequencies can be available during audiometry or the range can be abbreviated to eight frequencies. The default setting is **Aud Range Normal**. To select the abbreviated frequency range:

Position the square cursor in front of the feature **Aud Range Narrow**. Press the  button to save this narrow range for audiometric testing.

The word **SAVED** will appear in the lower right-hand corner and the asterisk now appears in front of the narrow range selection. The normal range of frequencies includes 125 Hz through 8000 Hz. The narrow range of frequencies includes 500 Hz through 6000 Hz. In the **AUD**


mode, if the narrow range is selected, the  and  **Hz** buttons will allow you to scroll through this abbreviated frequency range only. Both the screen and printout will still be labelled with the full range of frequencies (i.e., 125 Hz through 8000 Hz).

## PRINT – AUDIOGRAM / PRINT – AUD TABLE

The audiometric test results can be printed out in an audiogram format (**PRINT - AUDIOGRAM**) or in a tabular format (**PRINT - AUD TABLE**). The default setting for this function is the print audiogram format.

**NOTE:**When a specific frequency is not tested, the result will be a break in the audiogram on the printout. This eliminates the assumption that a threshold exists at that untested frequency.

To change the print option, move the cursor in front of the description PRINT - AUD TABLE.

Press the  icon to save this format as the new default parameter. The word **SAVED** appears in the lower right-hand corner of the display to indicate that this new setting has been saved.

With the **PRINT - AUD TABLE** selected, all audiometric test results will appear in a table with the frequency range typed horizontally along the top of the table followed by two lines of test data. The test results for the right ear will appear next to the letter **R** and below each frequency tested. Similarly, the test results from the left ear will follow below the right ear results.

**NOTE:**The **PRINT - AUD TABLE** setting selects the format for the printout only. An audiogram always appears on the LCD while in **AUD Mode**.

## DEF XDUCER TDH39 / DEF XDUCER INSERT

The **TDH39 HEADPHONES** are the factory default transducers. To select the **INSERT EARPHONES** as the default start up option, move the cursor next to the **DEF XDUCER**

**INSERT** selection and press  to save. **SAVED** will appear on the bottom right corner of the LCD. An **“\*”** will appear next to the **DEF XDUCER INSERT** option to denote the selection.

## Program Menu Page 2 option descriptions

### DATA XFER CONFIG

* 115.2 KBAUD	* NO PARITY + 8-BIT
57.6 KBAUD	ODD PARITY + 7-BIT
38.4 KBAUD	EVN PARITY + 7-BIT
17.2	SPC PARITY + 7-BIT
9600 BAUD	* XON/XOFF DISABLED
4800 BAUD	XON/XOFF ENABLED

These settings are used to allow the data transfer from the GSI 39 to a computer. The settings must be the same on the GSI 39 and the computer. The factory defaults are defined by an \*.

### POWER UP SETTINGS

TYMP  
\* TYMP REFLEX  
AUDIO

The feature allows you to program the GSI 39 to start up into any of the three test modes.



An asterisk will denote which option is selected to appear on the LCD when the system is first powered up. The factory default setting is TYMP REFLEX.

---

## PRN HEADER GSI/PRN HEADER OFF/PRN HEADER CUSTOM

There are three options to choose how the print header will be handled:

### PRINT HEADER GSI

This is the factory default setting for this feature. Each time the **Print Screen**  or **Print All**  tests in memory buttons are pressed, the printout will begin with the label **GSI 39**.


### PRINT HEADER OFF

If this option is selected, no header will be printed before any test results, which will save space and printout time.



### PRINT HEADER CUSTOM


Select this option to design a custom header, which might be the name of your own facility, department or company.

To type in the custom header, position the square cursor in front of **PRN HEADER**

**CUSTOM**. Press  to select it as the new default setting. The word **SAVED** appears in the lower right corner.

If **PRN HEADER CUSTOM** is selected, a line cursor will flash in the left-hand corner below the words **PRN HEADER CUSTOM**. To “type” in the desired header, use the dB HL knob. Rotating the knob clockwise will sequence you through the alphabet in the forward direction and rotating this knob counterclockwise will sequence you through the characters in a reverse direction. The available character set is: A -Z; 0 - 9; and a blank space. A total of 35 character spaces are available. Once you have cycled through to the


desired character, press the  button to store it. The cursor will move into position for the next character. Select the next character and press  to store. When the custom header is complete, press the **PROG** button to exit from the submenu.

To change/delete a previously saved character, press the  to position the cursor at that character. Use the HL knob to select the new character to change or select the blank space to delete.

**NOTE:**To center the header, consider the length of the name to be inserted and calculate from the left margin where you would like the header to begin. Type blank spaces to the start point of your custom header. If you begin to enter the characters for the header from the left margin, the header will be printed from the left margin on the printout.

## **INTERNAL PRINTER/EXTERNAL PRINTER**

These items toggle between either printing to the internal printer (4" paper) or sending the information to an external printer. The external printer is connected through a USB port on the back panel. The printer must be a DeskJet with PCL3 or PCL3GUI protocol.

To select the printer, move the cursor next to the Internal or External printer and press the  button to save the setting.



## RESET TO DEFAULTS

This option will reset the programmable settings to the Cardinal Health/GSI factory defaults.

Program Mode - Factory Default Settings	
AUDIOMETRY RESULTS	- PRINT - AUDIOGRAM
REFLEX RESULTS	- REFLEX HL + CURVE
226 Hz Reflex	- Ipsi 1000 Hz
1 kHz Reflex	- Ipsi 2000 Hz
NORMAL BOX	- NORMAL BOX ASHA
NEWBORN NORM	- NEWBORN NRM ON 1k
50th PERCENTILE	- 50th PERCENT ON 1k
BASELINE	- BASELINE OFF 1k
AUTOSTART	- AUTOSTART OFF 1k
AUDIOMETRY RANGE	- AUD RANGE NORMAL
DEFAULT TRANSDUCER	- TDH 39
PRINT HEADER	- PRN HEADER GSI
LANGUAGE	- ENGLISH
DATA XFER CONFIG	- 115.2 KBAUD
	NO PARITY + 8-BIT
	XON/XOFF DISABLED
POWER UP SETTING	- TYMP REFLEX
	226 Hz
PRINTER TYPE	- INTERNAL PRINTER
AUTO HL SETUP	
TEST FREQ (Hz)	- 500 Hz
	- 1000 Hz
	- 2000 Hz
	- 3000 Hz
	- 4000 Hz
	- 6000 Hz
INTEN RANGE (DBHL)	
MIN DB	- 0
MAX DB	- 90
START TEST EAR	- RIGHT
SCORING RULE	- 2 OUT OF 3
TONE FORMAT	- STEADY

### **Exiting the program mode**

Press the **PROG** button to exit the program mode and return to the previously selected test mode.

# Chapter 5


# Routine Maintenance

---

Blank page.

## **PreTest Tymp checks**

For your convenience, a test cavity is provided with your instrument. This test cavity enables you to quickly verify, on a daily basis, the proper calibration of your unit. Cardinal Health/GSI strongly recommends that you make this quick check a part of your daily routine.



**Figure 1: Test cavity.**


## Calibration Quick Check for 226 Hz Probe

To initiate the quick check, select the Tymp only mode and insert the probe into the 0.5 cm<sup>3</sup> opening on the test cavity. See Figure 1.

The instrument is designed to start automatically, it is important that the probe is inserted as quickly and as smoothly as possible. During the calibration check, the probe must be held carefully and without movement. Do not place the probe on the same counter as the instrument or any moving object during this check as mechanical noise may be picked up by the probe and interfere with the calibration check.


The calibration check will start automatically if the probe has been inserted into the cavity properly. This is confirmed by the green lamp changing from blinking to a steady condition. If the **orange** lamp is illuminated, the probe is not properly positioned within the cavity so that a large pressure leak exists. If the **yellow** lamp is illuminated, the probe tip has been occluded. In either case, remove the probe and wait for the blinking **green** lamp. Insert the probe once again. If necessary, clean the probe tip as described later in this chapter.

The green lamp will resume blinking when the probe is removed from the test cavity. The tympanogram on the display represents the response from the 0.5 cm<sup>3</sup> hard walled cavity. The ECV (ear canal volume) should read 0.5. The letters NP will appear alongside the pressure (daPa) and compliance (cm). Three dashed lines - - - will appear alongside the gradient (GR)

Using the same sequence, place the probe in the test cavity opening labelled 2.0 cm<sup>3</sup>. The resulting tympanogram should be identical other than the ECV should read 2.0 cm<sup>3</sup>. The same sequence can be followed with the 5.0cm<sup>3</sup> opening on the test cavity. To keep a record of this test cavity calibration check, simply press the  button on the front panel of the instrument.

Since sound pressure will vary with altitude and barometric pressure, some variation from the 0.5, 2.0 and 5.0 cm<sup>3</sup> readings may be observed. Your instrument is carefully calibrated at our factory, which is at approximately 850 feet above sea level. If you are located at an elevation of 1000 feet or higher, your instrument may need to be recalibrated to account for your elevation (see *Altitude Adjustment* in this chapter for more details). It is not necessary to recalibrate for barometric pressure changes on a daily basis. Keep in mind that a change in barometric pressure (i.e., from low to high or vice-versa) will slightly affect the test cavity readings.

### Calibration Quick Check for Combo Probe

To perform a calibration quick check with the 226 Hz probe tone using the Combo probe, follow the directions on the previous page. To initiate the quick check for the 1000 Hz probe tone, select the Tymp only mode and 1000 Hz probe tone from the front panel. Insert the probe into the 0.5 cm<sup>3</sup> opening on the test cavity. See Figure 1 earlier in this chapter. If the Autostart option is **Off**, press the  to begin the measurement

**NOTE:**The factory default settings on the Autostart option of the 1000 Hz probe tone will be set to **Off**.

When the test begins, the green lamp should change from blinking to a steady condition. If the blinking orange lamp is illuminated, the probe is not properly positioned within the cavity so that a large pressure leak exists. If the steady orange light is illuminated after the tip is inserted securely into the cavity, there may be an occlusion as the system considers the measurement too small to begin the test. In either case, remove the probe and wait for the blinking green lamp. Insert the probe once again. If necessary, clean the probe tip as described later in this chapter.

When the test sequence is completed, the green lamp on the probe is no longer illuminated. The green lamp will resume blinking when the probe is removed from the test cavity. If baseline is **Off** (factory default setting), you will see the flat line across the display at the amplitude of the cavity reading. If baseline is **On**, then you should see a flat line along the bottom of the display at the 0 value. The C1 value for 1000 Hz probe tone is represented in mmhos and not converted to a volume reading. The reading for the 0.5 Cavity should be 2.2 mmhos.





Figure 2 - 1000 Hz probe tone in 0.5 cavity with the baseline turned off.

## **GSI 39**

---

Since the test cavity is a hard-walled cavity, the tympanogram should be a flat line indicating that there is no mobility in the system. The instrument places the letters NP next to the mmho and daPa headers to indicate that there is no peak compliance and, therefore, no peak pressure can be determined during the quick check.

Using the same sequence, place the probe in the test cavity opening labeled 2.0 cm<sup>3</sup>. Note that when the measurement is finished, the display will change to 10 mmho scale and the normative data will no longer be displayed. The C1 value for the 2.0 cm<sup>3</sup> cavity should be approximately 8.85 mmhos.

To keep a record of this test cavity calibration check, press the Print All  button on the front panel of the instrument.



### Altitude adjustment







The altitude calibration adjustment allows you to ‘correct’ the ear canal volume (ECV) measurement and test cavity volume measurement for variations due to altitude. The instrument is a pressure sensitive device that makes measurements relative to ambient air pressure. Changes in air pressure due to weather or altitude will affect the ECV readout of the instrument. The slight pressure change resulting from changing weather conditions will usually yield volume readouts with  $\pm 0.1 \text{ cm}^3$  of the expected cavity value, but pressure changes due to altitude can shift these cavity values by as much as 30%. These changes in pressure do not affect the accuracy of the compliance measurement system in any way. However, you may prefer that the instrument gives ECV values as they would appear at sea level. The altitude calibration mode allows you to adjust the Auto Tymp without the services of a qualified Cardinal Health/GSI representative.






226 Hz Probe		1000 Hz Probe	
Altitude in Feet	Equivalent 2.0 cc Reading	Altitude in Feet	mmho Reading
0	2.0 $\pm 0.1$	0	8.85 mmho $\pm 0.44$
1000	2.1 $\pm 0.1$	1000	9.20 mmho $\pm 0.46$
2000	2.2 $\pm 0.1$	2000	9.56 mmho $\pm 0.48$
3000	2.2 $\pm 0.1$	3000	9.91 mmho $\pm 0.50$
4000	2.3 $\pm 0.1$	4000	10.3 mmho $\pm 0.52$
5000	2.4 $\pm 0.1$	5000	10.6 mmho $\pm 0.53$
6000	2.5 $\pm 0.1$	6000	11.1 mmho $\pm 0.56$
7000	2.6 $\pm 0.1$	7000	11.5 mmho $\pm 0.58$
8000	2.7 $\pm 0.1$	8000	12 mmho $\pm 0.60$
9000	2.8 $\pm 0.1$	9000	12.4 mmho $\pm 0.62$
10000	2.9 $\pm 0.1$	10000	12.8 mmho $\pm 0.64$


Table 1: 226 Hz Probe Tone Altitude Correction.

Table 2: 1000 Hz Probe Tone Altitude Correction.

**NOTE:** Operation of the system at an altitude of 10,000 feet may affect the ability to pressurize to the maximum 5.00  $\text{cm}^3$ .

To enter altitude calibration, press ,  and  simultaneously. The LCD will now display the Settings main menu. To enter the Altitude mode, move the square cursor in position using  and  just to the left of the Altitude mode and press the PAGE  button.

1. When entering the altitude mode, the display will read as follows:  
Altitude Mode  
ECV 2.0  
cm<sup>3</sup> 9.99  
Standard
2. Select **226 Hz** or **1000 Hz** probe tone.
3. Place the probe into the 2.0 cm<sup>3</sup> cavity provided with the instrument and check the cm<sup>3</sup> value against the altitude correction table for accuracy.
4. If the measured volume is not within the published table value  $\pm 0.1$  cm<sup>3</sup>, then you should exit the altitude mode by pressing the **PROGRAM MODE** button and contact field service. Providing the measured volume agrees with the published table  $\pm 0.1$  cm<sup>3</sup>, you may proceed with the altitude adjustment.
5. With the probe still in the 2.0 cm<sup>3</sup> cavity, press the **PROG** button to enter the custom calibration mode. **CUSTOM** will appear on the fourth line of the display.
6. The value now displayed in the cm<sup>3</sup> display area is the volume measured and adjusted to the current altitude. If the value displayed is 2.0 cm<sup>3</sup>, the volume is adjusted to the current site. If the value is not 2.0 cm<sup>3</sup>.
7. press the  **SAVE** button to customize the volume measurement to the current altitude. The measured volume should now read 2.0 cm<sup>3</sup>.
8. To exit the altitude mode, press the **PAGE**  button to return to the Settings main menu.
9. Move the cursor using  and  to **Back to Normal** and press the **PAGE**  button to return to the Normal mode.

 **WARNING** Cardinal Health/GSI recommends only trained personnel enter the Calibration and Diagnostic submenus listed below the Altitude Adjustment in Calibration mode.

## **Pre-Test Audiometric Checks (Version 3 and 4 only)**

### **Noise recovery period**

Exposure to high levels of sound (e.g., unmuffled lawn mowers, loud music, gunfire) tends to create a temporary threshold shift (TTS) which diminishes with time after exposure. Any subject/patient tested soon after such exposure may exhibit a hearing loss that does not reflect his or her normal hearing threshold. It is, therefore, important that the testing procedure prescribe some time interval - usually at least 16 hours- between the last exposure to high-level sounds and the administration of any hearing test.

---

**Elimination of ambient noise**

Excessive noise in the test environment during audiometric testing such as that produced by conversation, typewriters, public address systems reduces test validity as it tends to mask the test signals, particularly at the lower frequencies where earphone cushions provide less effective attenuation. An acoustically treated room may be required if ambient noise reaches objectionable levels (i.e., sufficient to cause apparent hearing loss at the low frequencies). Also, Audiocups are available from Cardinal Health/GSI as an optional accessory. If the person being tested is in the same room as the audiometer, it is recommended that he/she be seated about three feet (1 meter) away from the instrument.

Maximum permissible noise levels are specified by the American National Standards - *Criteria for Permissible Background Noise during Audiometric Testing*, ears covered with earphones (S3.1 1991 revised). Table 3 shows the maximum background levels that can be present inside the room while a valid hearing test is being conducted. For more comprehensive information about hearing testing and hearing conservation, the user is referred to the Bibliography.



<b>Frequency (Hz)</b>	<b>Test Room Maximum dB SPL* in 1/3 Octave Band</b>
<b>125</b>	<b>29.0</b>
<b>250</b>	<b>17.5</b>
<b>500</b>	<b>14.5</b>
<b>750</b>	<b>16.5</b>
<b>1000</b>	<b>21.5</b>
<b>1500</b>	<b>21.5</b>
<b>2000</b>	<b>23.0</b>
<b>3000</b>	<b>28.5</b>
<b>4000</b>	<b>29.5</b>
<b>6000</b>	<b>33.0</b>
<b>8000</b>	<b>38.5</b>



**Table 3**

## Biological Check

For Tympanometry and Reflex tests, the best way to determine that your instrument is operating properly is to perform a daily check on a normal ear - your own ear if possible. This allows you to listen for the probe tone and the stimulus tone (during reflex) and to determine if the air pressure system is working properly. Keep a copy of your chart for a day-to-day reference in checking your instrument.

To perform a biological check in Audiometry, select the Audiometry (**AUD**) mode button. The display changes from the tympanogram format to an Audiogram format. Select **Headphone** or **Insert Phone**. (When changing transducers, the icon for the new transducer will flash on the

LCD until the button is pushed again.) The  and  **Hz** buttons allow you to select each frequency and the **dB HL** knob allows you to alter the intensity of each frequency. Position the test headset on your head so that each earphone is covering the appropriate ear (i.e., **red** is right and **blue** is left). Select the right earphone by pressing the front panel button labelled **R** and check for the following while depressing the **Present bar**:

- a. Pressing the  **Hz** button causes the frequency to change to a lower frequency.  
Pressing the  **Hz** button causes the frequency to change to a higher frequency.
- b. Each frequency or tone is pure (i.e., there is no distortion or crackling sound present).
- c. Rotating the **dB HL** knob in a clockwise direction causes the tone to increase (become louder) in intensity. Rotating the **dB HL** knob in a counterclockwise direction causes the tone to decrease in intensity (become quieter).

Since individual thresholds can shift up or down as much as 5 dB from one day to the next, variation within this range may be considered acceptable. Variations that exceed this range, however, are likely to reveal problems that require attention. The routine maintenance checks described in this chapter, may suggest the source and solution to the problem. If they do not, the instrument should receive technical service by a Cardinal Health/GSI certified technician before further use.

## **Preventive Maintenance**

Preventive maintenance does not require access to the interior of the instrument and may be performed by the you.

For the GSI 39, preventive maintenance consists of periodically cleaning and inspecting the exterior of the instrument. We also recommend cleaning and inspecting the accessories such as the probe and/or earphones. It is recommended that you develop a schedule for these purposes.

### **Cleaning the system**

Turn **OFF** the system power before cleaning the instrument. Do not permit solutions or sterilization agents to seep into the electronic portions of the system. Take special care around controls, connectors and panel edges. Do not use any abrasive cleaners.

Remove any dust from the exterior of the system with a soft brush or cloth. Use a brush to dislodge any dirt on or around the connectors and panel edges. Remove stubborn dirt with a soft cloth slightly dampened with a mild detergent solution or cold sterilization agent.

### **Recommended cleaning solutions**

Housing components should be wiped with a damp cloth containing soap and water, isopropyl alcohol, ammonia based cleaners or bleach based cleaners.

### **Cleaning patient contact reusable devices**

To help ensure patient safety, prevent cross infection and provide effective service, Cardinal Health/GSI patient contact devices must be properly maintained. Maintenance should include cleaning before each use.

The earphone cushions and patient hand switch can be wiped with a slightly damp cloth containing soap and water, isopropyl alcohol, ammonia based cleaners or bleach based cleaners. Gently wipe the earphone cushions with the slightly damp cloth taking care not to get moisture in the speaker portion of the earphones.

We recommend eartips be discarded after use. We do not recommend cleaning and re-using the rubber probe eartips for tympanometry or the foam eartips used with the insert earphones.



**WARNING** It is recommended that all repairs be performed by a qualified Cardinal Health/GSI service representative only. You have the sole responsibility for any malfunctions resulting from improper maintenance or repair by anyone other than an authorized Cardinal Health/GSI representative.


## Probe care - 226 Hz Probe

With use, cerumen can work its way up inside the probe nose cone (probe tip). During the warm-up period each day and throughout the day, inspect the probe tip to make sure it is clean and free of cerumen. Refer to the following instructions for cleaning and maintaining the instrument's probe.

### Probe nose cone cleaning

Remove the nose cone portion of the probe:

1. Hold the body of the probe in one hand (e.g., left) near the tip and grasp the nose cone of the probe in the other hand (e.g., right).
2. Rotate the nose cone portion of the probe counterclockwise until the nose cone is completely separated from the probe (Figure 3).
3. Place the probe body securely on a table and inspect the nose cone for cerumen. Use a pipe cleaner to remove any cerumen by inserting the pipe cleaner through the back portion of the nose cone and pulling it through the front opening. It may be necessary to repeat this several times to remove all the cerumen.




**Figure 3: Probe nose cone removal.**

**NOTE:** The probe nose cone can be sterilized via many conventional methods including autoclaving.



### The O-Ring

There is an O-Ring seated at the end of the threads on the probe. As a preventative maintenance measure, and to ensure that the nose cone of the probe unscrews easily, do not clean or remove the lubricant from the O-Ring. If the O-Ring appears to be void of any lubricant, or if the nose cone itself was difficult to remove, apply a high-quality synthetic lubricant such as those considered “food-grade.” Refer to Figure 4 and apply as described in the instructions that follow.



**Figure 4: O-Ring care.**

**A:** Cotton swab.

**B:** Lubricant

**C:** O-Ring (enlarged for detail).

1. Place a small drop of lubricant at the front outer surface of the O-Ring.
2. Using your finger or a cotton swab, spread a thin layer of lubricant completely around the front and outer surface of the O-Ring. Ensure that no lubricant spreads into the threaded area of the nose cone. Only a thin layer of lubricant is necessary. Excessive application or build-up may affect test results.

## The probe wire

Inside the probe body there is a metal tube that contains a wire required for cleaning purposes.

1. Carefully remove this wire from the metal tube (Figure 5). This will pull any cerumen out of the metal tube.



**Figure 5: Probe wire removal.**

2. Examine the wire for cerumen.
3. If necessary, clean the wire with a lint-free tissue.
4. Reinsert the wire into the metal tube and push it in as far as it can go.

**NOTE:**The wire must be inserted into the metal tube for the instrument to function properly.

## Probe reassembly

After cleaning, reassemble the probe nose cone to the probe body by screwing the cone back onto the probe. Take care to align the threads on both the probe body and the nose cone before screwing the pieces together. Only screw the nose cone on until it is finger tight. You might find it helpful to gently squeeze the two sides of the probe case together while screwing the nose cone into place.


**NOTE:**The probe nose cone must be screwed firmly in place to guard against any air leaks.

## **Probe Care - Combo Probe Tip**

To ensure measurement accuracy, it is essential to clean the probe tip daily to be certain that the tubes are clean and free of cerumen. If the lights on the probe are indicating an occlusion, cleaning the probe tubes will most likely rectify the situation.

The Cleaning Floss Kit (2000-9610) contains 2 sizes of floss, which can be used to clean the three metal tubes on the probe tip.

1. Remove the probe eartip and tygon tubing attached to the three metal probe tubes at the rear of the probe tip.



### **WARNING**

**Do not alter the length of the tubing by cutting it. When reconnecting the tygon tubing to the probe tubes, ensure that there are no sharp edges or burrs on the probe tubes that could cut the tygon tubing.**


## **GS1 39**

---

2. Using the smaller sized floss, for one of the smaller size metal tubes, insert the floss into the base of the black probe tip and pull the floss through the metal tube. Discard the used floss.
3. Repeat with a second piece of the small floss for the other small tube.
4. Use the larger floss to clean the large tube of the probe tip in the same manner. Do not reuse the cleaning floss.



5. Reconnect the tygon tubing to the metal probe tubes. The center tube has the larger diameter.



Tygon tubing should be replaced if debris can be seen in the tubing as this could affect accuracy of the measurements.

### **WARNING**

**Avoid getting the probe moist. Do not use the probe tip if it is wet or damp because the moisture may make its way to the sensitive electronic equipment at the end of the tygon tubing.**

## **Earphone Care (Versions 3 and 4 only)**

With proper care, the earphone and cords provided with the instrument (versions 3 and 4) should last a long time. Moisture should not be allowed anywhere near the earphone itself as this will damage the diaphragm and grill cloth, requiring its replacement. The earphone cushions can be wiped with a slightly damp cloth, taking care not to get moisture in the speaker portion of the earphones.

With extended use, earphone cords tend to fray internally at the connectors (i.e., between the cord and the instrument's connector, and between the cord and the earphone connector). This fraying may cause a decrease in the signal level or cause the signal to be intermittent. To check for this:

1. Position the test headset over your ears and select a frequency (e.g., 1000 Hz) at 35 dB HL.
2. Select the right earphone and press the **Present bar**.
3. While the Present bar is depressed, flex the earphone cord next to the connector at both ends.
4. Listen for an intermittent signal, an abrupt change in signal intensity level or a scratchy sound superimposed over the selected frequency that coincides with the flexing of the cord. The presence of any of these conditions indicates that the cord should be replaced.
5. Also, examine the earphone cord for cuts or tears in the covering shield and the earphone cushion for signs of damage. If either problem is noticed, the earphone cord or cushion should be replaced. Both parts are easily replaced without the need for recalibration. However, if the earphone receives shock damage or is replaced for any reason, the instrument will need to be recalibrated.
6. Repeat this same sequence with the left earphone.

### **Paper supply**

To streamline each testing session, it is a good idea to check the amount of paper left inside the printer compartment. Extra rolls of paper should be kept nearby.

**NOTE:** The number of tests per roll of paper will vary with the version Auto Tymper being used and the type of tests being performed. See *Printer Description* in the *Specifications* section of this guide for approximations. Replacement paper can be purchased from your local Cardinal Health/GSI Distributor or from the factory.

# Chapter 6

# Test Results

---

Blank page.



## Ear Canal Volume - 226Hz Probe Tone

### Normal

As a general rule, values for ear canal volume should be between 0.2 and 2.0 cm<sup>3</sup>. However, the normal values will vary with age and bone structure.

### Abnormal

An ear canal value of less than 0.2 cm<sup>3</sup> indicates an abnormal condition. If the probe is partially plugged with cerumen or if the probe is positioned up against the ear canal wall, a smaller than expected value will be measured. Also, if an individual has a relatively large bone structure for his/her age group and a smaller than expected value is measured, the probe could also be partially occluded or up against the canal wall. It is also possible to collapse the canal if the probe is held too firmly against it. Examine the Tympanogram and the reflex results to confirm your suspicions. If they are abnormal as well, it is good practice to repeat the test.

An ear canal volume greater than 2.0 cm<sup>3</sup> also may indicate an abnormal condition. An important application of the ear canal volume measurement is to determine if there is a perforation of the tympanic membrane. If there is a perforation due to trauma or due to the presence of a pressure-equalization (P-E) tube, the measured ear canal volume will be much larger than normal since the combined volume of the ear canal and the middle-ear space is being measured. The maximum ECV is 5.0 cm<sup>3</sup>, any space larger than that will be recorded as 5.0 cm<sup>3</sup> or may not seal.

## Compliance Peak

### Normal

The range of normals for compliance is 0.2 cm<sup>3</sup> to approximately 1.4 cm<sup>3</sup>. Some groups use a larger range up to 1.8 cm<sup>3</sup>. A measured compliance peak within this range indicates normal mobility within the middle-ear system.

### Abnormal

A compliance value of less than 0.2 cm<sup>3</sup> indicates a pathological condition as the middle-ear system is stiffer than normal. To distinguish the probable cause of the stiffening, the pressure value where this stiffened compliance peak occurs needs to be considered. For example, normal pressure along with a stiff middle ear system is indicative of otosclerosis, a severely scarred tympanic membrane or a layer of plaque across the tympanic membrane. On the other hand, abnormal pressure along with a stiffened middle-ear system is consistent with a poorly functioning eustachian tube with possible effusion (serous otitis media) or “glue ear.”

**NOTE:** If the measured compliance value is less than 0.1 cm<sup>3</sup>, the letters NP will be printed next to the heading cm<sup>3</sup> on the screen and printout. The letters NP indicate a poorly defined or flat Tympanogram. The Tympanogram may depict a very shallow peak.

A compliance value greater than 1.4 cm<sup>3</sup> (or 1.8 cm<sup>3</sup>) indicates a hyperflaccid tympanic membrane or a possible disarticulation depending upon how far above the normal range the value is. Generally speaking, a compliance value of greater than 3.0 cm<sup>3</sup> is indicative of a disarticulated ossicular chain. Further testing is necessary to confirm this suspicion.

**NOTE:** If a compliance value is measured to be greater than 1.5 cm<sup>3</sup>, the instrument automatically changes the range assigned to the graph to 3.0 cm<sup>3</sup>.

The validity of tympanometry and acoustic reflex testing is dependent upon a healthy tympanic membrane. A pathological condition at this membrane can mask the true condition of the middle ear.

## Pressure Peak

### Normal

Strict rules for middle-ear pressure indicate a normal range of  $\pm 50$  daPa. However, for most applications, a normal range of -150 daPa to +100 daPa is used.

### Abnormal

Very rarely will you obtain an extreme positive pressure condition. Some researchers have reported high positive pressures at the onset of acute otitis media.

Pressure values more negative than -150 daPa are indicative of a poorly functioning eustachian tube. The severity of this condition is determined by how negative the pressure is and its impact on the compliance peak.

If no pressure peak is measured over the pressure range of +200 daPa to -400 daPa, then the letters NP will appear on the screen and the printout. This indicates that no pressure peak was detected over this pressure range.

## Gradient

### Normal

When testing a child, the normal range for the gradient is between 60 and 150 daPa. (Infants may show higher gradient values due to the mobility of their ear canals.) The range of normal is somewhat narrower for adults (i.e., 50 to 110 daPa).

### Abnormal

A high gradient value (greater than the high end of the normal range per age group) is indicative of middle-ear effusion. The reduced compliance values and negative middle-ear pressure characteristic of developing or resolving otitis media with effusion (OME) will be manifested in a broad tympanogram with a large gradient value. However, abnormal gradient values may also be found in the absence of abnormal parameters. This could indicate a transient OME, so a retest after several weeks may be recommended.

When the middle ear's mobility is reduced to near 0 cm<sup>3</sup>, due to viscous effusion or a "glue-ear" condition, no gradient value can be measured. In this case, dashes (- -) will be displayed next to the letters GR.

Very low gradient values are associated with a flaccid middle ear system. These low values should be taken into consideration with the ear canal volume and compliance peak values to determine the probable use of the flaccid condition.

## Acoustic reflex

### Normal

For screening purposes, an ipsilateral or contralateral reflex measured at any one of the levels available per frequency can be considered normal. Obviously, the lowest values are desired. However, without knowing the hearing threshold level of the individual per frequency, it is difficult to make a more definite statement. Generally speaking, the reflex is reported to occur at between 70 and 90 dB HL above the hearing threshold in normals. Remember that these values apply to reflex threshold measurements and that your instrument does not permit reflex threshold measurements due to the use of a hand-held probe. The presence of a reflex in the absence of a compliance peak suggests that the tympanometric results should be considered invalid and the test repeated. This is true because if there is no compliance measured during tympanometry, it is not possible to measure any stiffening affect during the reflex stimulus presentation.

### Abnormal

If a pressure leak occurs during the reflex testing and the pressure system is unable to correct for this leak, the reflex test sequence is aborted. When this occurs, the test results are assigned the letters NT (Not Tested).

If no response is obtained at the third and final stimulus level, the instrument will indicate this with the letters NR or No. More detailed testing at the frequency where this occurred is required to determine the reason for the no response.

## Audiometry

### Normal

A normal response from a child should be at or below 20 dB HL. A normal response from an adult will be somewhat higher at or below 25 dB HL. Remember that these normal values assume a quiet environment during testing.


### Abnormal

In children, a failure to respond to a 20 dB HL (or lower) stimulus presentation during a retest performed four to six weeks after the initial test would indicate the need for more extensive diagnostic testing to determine the cause.

In adults, a failure to respond at or below 25 dB HL when the room noise levels are low indicates the need for more evaluation. However, the age and employment history of the individual must also be considered.

## Special Messages and Error Codes


Error code numbers and other special messages may be displayed on the screen or on the printout. These messages appear whenever an instrument error occurs or, in some instances, to apprise the operator of certain situations. For example, if there is no test result on the screen and

the **Print Screen**  button is pressed, the printer will indicate “No Test To Print”


Error codes will appear as a two-digit number prefixed by the letter “E”. If an error code appears, please repeat the operation that caused the error code to appear. If the error code appears for the second time, make a note of it and contact your Cardinal Health/GSI Service Representative, giving him/her the exact error code number.

**Sample Test Results**


Figures 1 through 8 illustrate test results from sample GSI 39 Auto Tymp printouts. The smoothness of the tympanogram tracing is determined by the amount of movement during the testing. Little or no movement during the testing provides a smoother tracing. Moving, talking or crying during testing leads to a more erratic looking tracing but does not dramatically affect the test results.




**Figure 1: 226 Hz Range of Normals**




**Figure 2: 226 Hz Abnormal Tymp**




**Figure 3: 226 Hz Abnormal Tymp**



**Figure 4: 226 Hz Abnormal Tymp**



**Figure 5: 226 Hz Abnormal Tymp**



**Figure 6: Infant Normative Data**

**Infant Normative Data Legend for Figure 6**

	<b>C1 y + 200 daPa</b>	<b>mmho y Peak Values</b>	<b>daPa Pressure at Peak</b>
<b>A</b>	5% 0.8	5% 1.2	5% -133
<b>B</b>	50% 1.4	50% 2.5	50% 0
<b>C</b>	95% 2.2	95% 4.8	95% 113




Figure 7: 1000 Hz Abnormal




Figure 8: 1000 Hz Abnormal



# Chapter 7


# Computer Interface

---


Blank page.

## Introduction

The Computer Interface provides the capability of transferring stored test results from the instrument to an external computer or data collection device via a USB connection.

 **WARNING** Accessory equipment connected to the analog and digital interfaces must be certified to the respective IEC standards (IEC 950 for data processing or IEC 60601-1 for medical equipment). Furthermore, all configurations shall comply with the system standard IEC 60601-1-1. Everyone who connects additional equipment to the signal input or signal output part configures a medical system, and is, therefore, responsible that the system complies with the requirements of the system standard IEC 60601-1-1. If in doubt, consult the technical service department or your local representative.


## Operation

Press the  button to transfer test results stored in memory. During data transfer, the message **DATA TRANSFER** will appear on the LCD screen.

### Transferring during normal operation


During normal testing operation, pressing the  button will transfer all stored test results sequentially.

### Transferring from memory pages

If the **PAGE** button is used to review individual test results stored in any of the 12 memory locations, the  button will transfer only the currently displayed stored test results. There is one exception to this rule: If the last (most recent) test result is displayed, the instrument assumes normal testing operation, and transfers all test results.

## Other LCD screen messages

### INVALID SELECTION


This message appears if the  button is pushed during any of the following circumstances:

- During presentation of an audiometric tone
- During a tympanometry test
- During a reflex test
- During Printing

### NO DATA AVAILABLE


This message appears if the  button is pressed and no results are stored.

### NOT AVAILABLE

This message appears if the  button is pressed and the computer is not properly connected.

## Data Transfer Program Mode

The Data Transfer Program mode is used to modify the GSI 39 USB interface configuration parameters to match the computer's USB Port settings.

Enter the Program mode by selecting the **PROG** button. Move the cursor to **Data Xfer Config** and press the  button to enter the submenu.

The following screen appears the first time the Data Transfer Program mode is entered showing the factory default settings.

```




PROGRAM MODE - DATA TRANSFER
* 115.2 kBAUD   * NO PARITY + 8-BIT
  57.6 kBAUD   ODD PARITY + 7-BIT
  38.4 kBAUD   EVN PARITY + 7-BIT
  19.2 kBAUD   SPC PARITY + 7-BIT
   9600 BAUD   * XON/XOFF DISABLED
   4800 BAUD   XON/XOFF ENABLED



```

These selections fall into three groups:

- Baud rate
- Parity and data bits
- Flow control

The default setting for each group has an asterisk (\*) before it so that it is easy to scan the settings for each group.

Selecting default settings for any of the groups is achieved in the same manner as the Program mode. Use the  or  buttons to move the solid square cursor down or up to the setting that you wish to select. Press the  button. The word **SAVED** will appear in the lower right corner of the screen and the asterisk (\*) will appear in front of the new setting.

To exit the Data Transfer Program submenu, move the cursor to the  press the  button. This will return you to the Program mode which can be exited by selecting **PROG**.

## **Computer Interface**

### **Interface configuration**

The configuration of the GSI 39 computer interface must be set to match the interface configuration of the computer. The GSI 39 defaults to 115 kBaud, no parity, 8 data bits, 2 stop bits and no communications flow control. The default settings for the baud rate, parity, number of data bits and flow control may be modified using the *Data Transfer Program Mode* explained earlier in this chapter.

### **Cable connections**

The computer interface provides a serial interface consisting of a USB connector.

Chapter 8  
**Electromagnetic Compatibility (EMC)**

---

Blank page.



## GSI 39 - Electromagnetic Compatibility (EMC) Information

Portable and mobile RF communications equipment can affect the GSI 39. Install and operate the GSI 39 according to the EMC information presented on this page and the next 4 pages.

The GSI 39 has been tested for EMC emissions and immunity as a standalone instrument. Do not use the GSI 39 adjacent to or stacked with other electronic equipment. If adjacent or stacked use is necessary, the user should verify normal operation in the configuration.

The use of accessories, transducers and cables other than those specified, with the exception of servicing parts sold by Cardinal Health as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the device.

Anyone connecting additional equipment is responsible for making sure the system complies with the IEC 60601-1-2 standard.

<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>		
<b>The GSI 39 is intended for use in the electromagnetic environment specified below. The customer or the user of the GSI 39 should assure that it is used in such an environment.</b>		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The GSI 39 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B Limits	The GSI 39 is suitable for use in all commercial, industrial, business, and residential environments.
Harmonic emissions IEC 61000-3-2	Class A Category	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

**GSI 39**

**Recommended separation distances between portable and mobile RF communications equipment and the GSI 39**

The GSI 39 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the GSI 39 can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the GSI 39 as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.23 \sqrt{P}$
<b>0.01</b>	0.12	0.12	0.23
<b>0.1</b>	0.37	0.37	0.74
<b>1</b>	1.17	1.17	2.33
<b>10</b>	3.70	3.70	7.37
<b>100</b>	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2** These guidelines may not apply to all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**GSI 39**

<b>Guidance and Manufacturer's Declaration - Electromagnetic Immunity</b>			
The GSI 39 is intended for use in the electromagnetic environment specified below. The customer or the user of the GSI 39 should assure that it is used in such an environment.			
<b>Immunity Test</b>	<b>IEC 60601 Test level</b>	<b>Compliance</b>	<b>Electromagnetic Environment-Guidance</b>
<b>Electrostatic Discharge (ESD)</b>  <b>IEC 61000-4-2</b>	$\pm 6$ kV contact  $\pm 8$ kV air	$\pm 6$ kV contact  $\pm 8$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be greater than 30%.
<b>Electrical fast transient/burst</b>  <b>IEC61000-4-4</b>	$\pm 2$ kV for power supply lines  $\pm 1$ kV for input/output lines	$\pm 2$ kV for power supply lines  $\pm 1$ kV for input/output lines	Mains power quality should be that of a typical commercial or residential environment.
<b>Surge</b>  <b>IEC 61000-4-5</b>	$\pm 1$ kV differential mode  $\pm 2$ kV common mode	$\pm 1$ kV differential mode  $\pm 2$ kV common mode	Mains power quality should be that of a typical commercial or residential environment.


# GSI 39

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The GSI 39 is intended for use in the electromagnetic environment specified below. The customer or the user of the GSI 39 should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance	Electromagnetic Environment-Guidance
<p><b>Voltage dips, short interruptions and voltage variations on power supply lines</b></p> <p><b>IEC 61000-4-11</b></p>	<p><b>&lt; 5% UT</b> (&gt;95% dip in UT) for 0.5 cycle</p> <p><b>40% UT</b> (60% dip in UT) for 5 cycles</p> <p><b>70% UT</b> (30% dip in UT) for 25 cycles</p> <p><b>&lt;5% UT</b> (&gt;95% dip in UT) for 5 sec</p>	<p><b>&lt; 5% UT</b> (&gt;95% dip in UT) for 0.5 cycle</p> <p><b>40% UT</b> (60% dip in UT) for 5 cycles</p> <p><b>70% UT</b> (30% dip in UT) for 25 cycles</p> <p><b>&lt;5% UT</b> (&gt;95% dip in UT) for 5 sec</p>	<p>Mains power quality should be that of a typical commercial or residential environment. If the user of the GSI 39 requires continued operation during power mains interruptions, it is recommended that the GSI 39 be powered from an uninterruptable power supply.</p>
<p><b>Power Frequency (50/60 Hz)</b></p> <p><b>IEC 61000-4-8</b></p>	<p>3 A/m</p>	<p>3 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

**Note:** *UT* is the a.c. mains voltage prior to application of the test level.

<b>Guidance and Manufacturer's Declaration - Electromagnetic Immunity</b>			
The GSI 39 is intended for use in the electromagnetic environment specified below. The customer or the user of the GSI 39 should assure that it is used in such an environment.			
<b>Immunity Test</b>	<b>IEC 60601 Test level</b>	<b>Compliance</b>	<b>Electromagnetic Environment-Guidance</b>
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the GSI 39, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.17 \sqrt{P}$ <p><math>d = 1.17 \sqrt{P}</math> 80 MHz to 800 MHz</p> $d = 2.23 \sqrt{P}$ 800 MHz to 2.5 GHz <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m)</p> <p>Field Strengths from fixed RF transmitters, as determined by an electromagnetic site survey, (a* on the next page) should be less than the compliance level in each frequency range. (b* on the next page)</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>
<p><b>Note 1:</b> At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p><b>Note 2:</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

## **GSI 39**

---

(a\*) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the GSI 39 is used exceeds the applicable RF compliance level above, the GSI 39 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the GSI 39.

(b\*) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

# Chapter 9

## Bibliography

---

American Speech-Language-Hearing Association (1990). "Guidelines for Screening for Hearing Impairment and Middle Ear Disorders". ASHA, 32 (Suppl.2), 17-24.

Criteria for Permissible Ambient Noise During Audiometric Testing (ANSI S3.1 - 1977).

de Jonge, R.R. (1986). "Normal Tympanometric Gradient: A Comparison of Three Methods". *Audiology*, 26, 299-308.

Koebshell, K.A. & Margolis, R.H. (1986). "Tympanometric Gradient Measured from Normal Pre-School Children", *Audiology*, 25, 149-157.

Margolis, R.H. & Heller, J.W. (1987). "Screening Tympanometry: Criteria for Medical Referral". *Audiology*, 26, 197-208.

Margolis, R.H. & Shanks, J.E., "Tympanometry". In Katz, J.(Ed.), *Handbook of Clinical Audiology*, Ed.3., Baltimore: Williams & Wilkins, 1985.

Margolis, R. H., Bass-Ringdahl, S., Hands, W., Holte, L. and Zapala, D. A. (2003) "Tympanometry in Newborn Infants - 1 kHz Norms" *Journal of the America Academy of Audiology*, 14 (7), 383-392.

Michael, P.L. and Bienvenue, G.R., "Noise Attenuation Characteristics of Supra-Aural Audiometric Headsets using the Models MX41/AR and 51 Earphone Cushions," *J.Acoust.Soc.Am.*, 70(5), Nov.1981, 1235-1238.

Methods for Manual Pure-Tone Threshold Audiometry (ANSI S3.21 1978). Newby, H.A., *AUDIOLOGY* (4th Ed.). New Jersey: Prentice-Hall Inc. (1979).

Paradise, J.L., Smith, C.G., Bluestone, C.D.(1976). "Tympanometric Detection of Middle Ear Effusion in Infants and Young Children", *Pediatrics*, 58 (2), 198-210.

U.S. Department of Labor, Occupational Noise Exposure, CFR 1910.95, March 8, 1983.

Blank page.



**Numerics**

- 1000 Hz REFLEX 4-10
- 226 Hz Probe Indicators 3-3
- 226 Hz REFLEX 4-10
- 50th PERCNT OFF 1k 4-6
- 50th PERCNT ON 1k 4-6

**A**

- Accessories supplied 2-3
- Acoustic Reflex 1-13
- Acoustic reflex 6-6
- Altitude adjustment 5-7
- Ambient noise 5-10
- atmospheric pressure 1-6
- Audiometric Checks (Version 3 and 4 only) 5-9
- Audiometric format during printing 3-20
- Audiometric Threshold 3-37
- Audiometry 6-6
- Audiometry Test Sequence (Versions 3 and 4 only) 3-31
- Audiometry testing information (Versions 3 and 4) 3-18
- Auto HL 1-13
- Auto HL mode 3-37
- AUTO HL SETUP 4-11
- Automated Audiometry (Auto HL) 1-13
- Automatic Hearing Level mode 3-37
- Automatic Hearing Level Procedure 3-39
- AUTOSTART OFF 1k 4-7
- AUTOSTART ON 1k 4-7

**B**

- BASELINE OFF 1k 4-6
- BASELINE ON 1k 4-6
- Bibliography 9-1
- Biological Check 5-11
- Bottom Panel 2-6
- Bottom Panel Label 2-6

**C**

- Cable connections 7-6
- Calibration Quick Check for 226 Hz Probe 5-4
- Calibration Quick Check for Combo Probe 5-5
- Catalog Listing 1-q
- Cautions 1-a
- Cleaning patient contact reusable devices 5-13
- Cleaning solutions 5-12
- Cleaning the system 5-12
- Combo Probe Indicators (226 Hz and 1000 Hz probe tone) 3-4
- Combo Probe Insertion 3-17
- Compliance Peak 1-13, 6-4

- Computer Interface 7-1, 7-6
- Connectors 2-5
- continuous bar 3-32
- contralateral 1-9
- Contralateral Acoustic Reflex 1-13
- Controls and Indicators 3-6
- Country Kits 1-r
- Customer responsibility 1-4

**D**

- daPa 1-6
- dashed bar 3-32
- Data Transfer Program Mode 7-5
- decaPascals (daPa) 1-6
- Display 2-4
- Disposal 1-12

**E**

- Ear Canal Volume 1-13
- Ear Canal Volume - 226 Hz Probe Tone 6-3
- ear canal volume (ECV) 3-22
- Earphone Care (Versions 3 and 4 only) 5-19
- ECV 3-22
- Elimination of ambient noise 5-10
- Error Codes 6-7
- European authority representative 1-e
- Exit audiometry 3-40
- Exit reflex 3-30
- Exit tympanometry/reflex 3-30
- Exiting the program mode 4-20

**F**

- FM 3-32
- Front Panel Controls and Indicators 3-6

**G**

- Glossary of terms 1-13
- GR 3-22
- Gradient 1-6, 1-8, 6-5
- gradient (GR) 3-22
- Green lamp 5-4
- GSI 39 Auto Tymp 1-3

**H**

- Helpful hints 3-13

**I**

- Indicators 3-6
- Initial set-up 2-8
- Inspection 2-3
- Instructing the patient/subject 3-18
- Interface configuration 7-6
- ipsilateral 1-9
- Ipsilateral Acoustic Reflex 1-13

## L

LCD screen messages 7-4  
Loading the paper 2-9

## M

Maintenance 5-12  
Manual mode 3-37  
Manual Threshold Audiometry 3-37  
Memory erase 3-41

## N

NEWBORN NRM OFF 1k 4-6  
NEWBORN NRM ON 1k 4-6  
Noise recovery period 5-9  
Normal Box 1-13  
NORMAL BOX ASHA 4-5  
NORMAL BOX OFF 4-5  
NT (Not Tested) 6-6

## O

Obtaining a seal 3-14  
OME 6-5  
Orange lamp 5-4  
O-Ring 5-15  
ossicular chain 1-7  
otitis media with effusion (OME) 6-5

## P

Panel Controls and Indicators 3-6  
Paper storage 2-10  
Paper supply 5-20  
Placement of earphones 3-19  
Placement of Insert earphones 3-19  
Preparing the probe assembly 3-4  
Pressure Peak 1-13, 6-5  
Pre-Test Audiometric Checks (Version 3 and 4 only) 5-9  
PreTest Tymp checks 5-3  
Printer and Display 2-4  
Printing history 1-d  
Printing test results 3-42  
Probe care - 226 Hz Probe 5-14  
Probe Care - Combo Probe Tip 5-17  
PROBE HZ 4-5  
Probe Indicators 3-3  
Probe nose cone cleaning 5-14  
Probe reassembly 5-16  
probe tip 3-17  
Probe Tone 1-13  
probe wire 5-16  
Program Mode 4-1, 4-3  
Programming the Auto HL Procedure 4-11

## R

Rear and Bottom Panel Labels and Connectors 2-5, 2-6

Rear Panel Labels and Connectors 2-5  
Recycling / disposal 1-12  
Reflex dB HL only 4-8  
Reflex dB HL plus curve 4-8  
REFLEX DISPLAY 4-8  
Reflex format 3-20  
Reflex Test Sequence 3-20  
Reflex yes/no 4-9  
Response handswitch (optional accessory) 3-19

## S

Safety notes 1-b  
Safety summary 1-a  
Sample Test Results 6-8  
Screening acoustic reflex 1-9  
screening acoustic reflex testing 1-10  
Screening Audiometry 1-13  
Screening audiometry 1-11, 3-36  
secretory otitis media 1-8  
Special Messages 6-7  
Specifications 8-1  
stapedial muscles 1-9  
stimulus 1-9

## T

Temporary programming of ipsilateral acoustic reflex test frequencies 3-26  
temporary threshold shift (TTS) 5-9  
Tests in memory 3-41  
Threshold Audiometry 3-37  
tip 3-17  
Transferring during normal operation 7-3  
Transferring from memory pages 7-3  
TTS 5-9  
TYMP OPTIONS 4-5  
Tympanogram 1-13  
Tympanometry and Gradient 1-6  
Tympanometry testing information 3-13  
Tympanometry/Reflex Test Sequence 3-20

## U

Unpacking and Inspection 2-3  
Upgrade 1-s  
Upgrade Kits 1-s

## W

Warnings 1-a  
Warranty 1-d  
WEEE - In Europe 1-12

## Y

Yellow lamp 5-4  
ympanometry 1-13